




PRODUCT CATALOG



 **PALTOP**

A KeystoneDentalGroup Brand

CONCEPTS

Internal Hex Connection

Narrow Ø 3/3.25 mm Workflow	6
Ø 3/3.25 mm Implants	7
Ø 3/3.25 mm Implants & Prosthetic Components	8
Ø 3/3.25 mm Prosthetic Components	9
Standard Ø 3.75/4.2/5 mm Workflow	20
Ø 3.75/4.2/5 mm Implants	21
Ø 3.75/4.2/5 mm Implants & Prosthetic Components	23
Ø 3.75/4.2/5 mm Prosthetic Components	24
Wide Ø 6 mm Workflow	36
Ø 6 mm Implants & Prosthetic Components	37
Ø 6 mm Prosthetic Components	38

Conical Connection

Conical Connection Workflow	48
Conical Connection Implants	49
Conical Connection Implants & Prosthetic Components	50
Conical Connection Prosthetic Components	51

LOCATOR®

LOCATOR® Workflow	62
LOCATOR® Prosthetic Components	63

Digital

Digital Abutment Level Workflow	66
Digital Implant Level Workflow	67
Digital Custom Milled Solutions Workflow	68
Digital Prosthetic Components	69

Surgical Kits & Additional Tooling

Premium Surgical Kit	76
Fully Guided Surgical Kit	82
Additional Guided Kit Tools	86
Basic Complete Surgical Kit	88
Mini Kits	92
Additional Tooling	100

Terms and Conditions


Terms and Conditions	102
----------------------	-----

Implant Specifications

Implant Specifications	104
------------------------	-----

Disclaimer: Some products may not be regulatory cleared/released for sale in all markets.

IH INTERNAL HEX CONNECTION




NARROW

PLATFORM

NARROW

WORKFLOW



Advanced Implant

ADVANCED CLASSIC	DESCRIPTION	REF #
	Ø 3.25 X 10 mm	20-70018
	Ø 3.25 X 11.5 mm	20-70019
	Ø 3.25 X 13 mm	20-70020
	Ø 3.25 X 16 mm	20-70021

Advanced + Implant

ADVANCED+ SLIGHTLY ACTIVE APEX	DESCRIPTION	REF #
	Ø 3 X 10 mm	20-70034P
	Ø 3 X 11.5 mm	20-70035P
	Ø 3 X 13 mm	20-70036P
	Ø 3 X 16 mm	20-70037P
	Ø 3.25 X 10 mm	20-70018P
	Ø 3.25 X 11.5 mm	20-70019P
	Ø 3.25 X 13 mm	20-70020P
	Ø 3.25 X 16 mm	20-70021P

Dynamic Implant

Dynamic	DESCRIPTION	REF #
	Ø 3 X 10 mm	21-70034
	Ø 3 X 11.5 mm	21-70035
	Ø 3 X 13 mm	21-70036
	Ø 3 X 16 mm	21-70037
	Ø 3.25 X 10 mm	21-70018
	Ø 3.25 X 11.5 mm	21-70019
	Ø 3.25 X 13 mm	21-70020
	Ø 3.25 X 16 mm	21-70021

PAI Implant

PAI	DESCRIPTION	REF #
	Ø 3.25 X 10 mm	29-70018
	Ø 3.25 X 11.5 mm	29-70019
	Ø 3.25 X 13 mm	29-70020
	Ø 3.25 X 16 mm	29-70021

PAI Total Coverage (TC) Implant

	DESCRIPTION	REF #
	Ø 3.25 X 10 mm	29-70018TC
	Ø 3.25 X 11.5 mm	29-70019TC
	Ø 3.25 X 13 mm	29-70020TC
	Ø 3.25 X 16 mm	29-70021TC

Cover Screw

	DESCRIPTION	REF #
	NP for Ø 3.25 mm Implant	80-70007

Concave Healing Caps

	DESCRIPTION	CUFF	REF #
	Ø 4 mm	1 mm	55-70009
	Ø 4 mm	2 mm	55-70010
	Ø 4 mm	3 mm	55-70011
	Ø 4 mm	4 mm	55-70014
	Ø 4 mm	5 mm	55-70015
	Ø 4 mm	6 mm	55-70018
	Ø 4 mm	7 mm	55-70019

Straight Healing Caps

	DESCRIPTION	CUFF	REF #
	Ø 4 mm	2 mm	55-70028
	Ø 4 mm	3 mm	55-70029
	Ø 4 mm	5 mm	55-70030

Impression Copings

	DESCRIPTION	REF #
	Open Tray Impression Coping	56-70016
	Open Tray Impression Coping, Slim	56-70033
	Closed Tray Impression Coping	56-70017
	Plastic Snap-On Closed Tray Impression Coping With Abutment	56-70035

Scan Abutment

	DESCRIPTION	REF #
	NP	30-70102

Implant Analogs

	DESCRIPTION	REF #
	Implant Analog	56-70015
	Digital Analog	30-70502

Titanium Temporary Abutments

	DESCRIPTION	CUFF	REF #
	Engage	0.5 mm	45-70031
	Engage	2 mm	45-70055
	Engage	3 mm	45-70057
	Non-Engage	0.5 mm	45-70032
	Non-Engage	2 mm	45-70056
	Non-Engage	3 mm	45-70058

Titanium Temporary Abutments, Immediate

	DESCRIPTION	CUFF	REF #
	NP	1.5 mm	45-70036
	NP	3 mm	45-70038
	Hex Driver for Ratchet for Immediate Temporary Abutment		60-70244

PEEK Concave Temporary Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	45-70020
	Straight	2 mm	45-70021
	Straight	3 mm	45-70022
	Angulated 15°	1 mm	45-70023
	Angulated 15°	2 mm	45-70024
	Angulated 15°	3 mm	45-70025

Fully Castable Straight Abutment

	DESCRIPTION	REF #
	Engage	46-70012
	Non-Engage	46-70017

Castable Abutments

	DESCRIPTION	REF #
	Gold Base Engage	46-70019
	Gold Base Non-Engage	46-70021
	Cobalt Chrome Base Engage	40-70078
	Cobalt Chrome Base Non-Engage	40-70178

Titanium Straight Abutments

	DESCRIPTION	CUFF	REF #
	Ø 4 mm	1 mm	40-70024
	Ø 4 mm	2 mm	40-70025
	Ø 4 mm	3 mm	40-70026
	Anatomic Ø 4 mm	1 mm	40-70053
	Anatomic Ø 4 mm	2 mm	40-70054
	Anatomic Ø 4 mm	3 mm	40-70055
	Anatomic Concave Ø 4 mm	1 mm	40-70056
	Anatomic Concave Ø 4 mm	2 mm	40-70057
	Anatomic Concave Ø 4 mm	3 mm	40-70058

Titanium Angulated Abutments

	DESCRIPTION	CUFF	REF #
	Anatomic 15°	1 mm	40-70027
	Anatomic 15°	2 mm	40-70028
	Anatomic 15°	3 mm	40-70029
	Anatomic Concave 15°	1 mm	40-70059
	Anatomic Concave 15°	2 mm	40-70060
	Anatomic Concave 15°	3 mm	40-70061
	Anatomic 20°	1 mm	40-70030
	Anatomic 20°	2 mm	40-70031
	Anatomic 20°	3 mm	40-70032

Titanium Angulated Abutments cont.

	DESCRIPTION	CUFF	REF #
	Anatomic Concave 20°	1 mm	40-70062
	Anatomic Concave 20°	2 mm	40-70063
	Anatomic Concave 20°	3 mm	40-70064
	Knife Edge 15°		40-70312

Snap-On Abutment System (SAS)

	DESCRIPTION	LENGTH	CUFF	REF #
	Ø 4.5 mm	4 mm	1 mm	44-70021
	Ø 4.5 mm	4 mm	2 mm	44-70022
	Ø 4.5 mm	4 mm	3 mm	44-70023
	Ø 4.5 mm	4 mm	4 mm	44-70024
	Ø 4.5 mm	4 mm	5 mm	44-70025
	Ø 4.5 mm	5.5 mm	1 mm	44-70026
	Ø 4.5 mm	5.5 mm	2 mm	44-70027
	Ø 4.5 mm	5.5 mm	3 mm	44-70028
	Ø 4.5 mm	5.5 mm	4 mm	44-70029
	Ø 4.5 mm	5.5 mm	5 mm	44-70030

Titanium Base Abutments Round

	DESCRIPTION	LENGTH	CUFF	REF #
	Engage	4 mm	0.5 mm	40-70215
	Engage	4 mm	1 mm	40-70322
	Engage	4 mm	2 mm	40-70295
	Engage	4 mm	3 mm	40-70217
	Engage	6 mm	0.5 mm	40-70219
	Engage	6 mm	1 mm	40-70324
	Engage	6 mm	2 mm	40-70297
	Engage	6 mm	3 mm	40-70221
	Non-Engage	4 mm	0.5 mm	40-70216
	Non-Engage	4 mm	1 mm	40-70323
	Non-Engage	4 mm	2 mm	40-70296
	Non-Engage	4 mm	3 mm	40-70218
	Non-Engage	6 mm	0.5 mm	40-70220
	Non-Engage	6 mm	1 mm	40-70325
	Non-Engage	6 mm	2 mm	40-70298
	Non-Engage	6 mm	3 mm	40-70222


Titanium Base Abutments Round Anterior

	DESCRIPTION	LENGTH	CUFF	REF #
	Engage	6 mm	0.5 mm	40-70302
	Engage	6 mm	2 mm	40-70303
	Engage	6 mm	3 mm	40-70304

Titanium Multi-unit Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	42-70011
	Straight	2 mm	42-70012
	Straight	3 mm	42-70013
	Straight	4 mm	42-70014
	Straight	5 mm	42-70015
	Angulated 17°	3 mm	42-70016
	Angulated 17°	4.5 mm	42-70017

Titanium Multi-unit Components

	DESCRIPTION	REF #
	Closed Tray Multi-unit Impression Coping	56-72003
	Open Tray Multi-unit Impression Coping	56-72002
	Plastic Snap-On Closed Tray Impression Coping	56-70049
	Plastic Snap-On Closed Tray Impression Coping, Ø4.5mm, (5 pack)	56-70062
 	Multi-unit Analog	56-72001
	Multi-unit Digital Analog	30-70504
 	Multi-unit Healing Cap, Tapered	55-72011
	Multi-unit Healing Cap	55-72010
 	Titanium Multi-unit Cylinder H 10 mm	42-70020
	Titanium Multi-unit Cylinder H 15 mm	42-70030
	Fully Castable Multi-unit Abutment	42-70021

Titanium Multi-unit Components cont.

	DESCRIPTION	REF #
	Titanium Retentive Multi-unit Interface Coping	42-70028
	Multi-unit Scan Abutment	30-70103

Titanium Multi-unit Screws

	DESCRIPTION	REF #
	Multi-unit Screw UNF 1-72	80-70005
	Open Tray Multi-unit Impression Coping Screw UNF 1-72	80-70028
	Abutment Screw UNF 1-72 Direct	80-70002*

Titanium Straight Single-unit Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	41-70151
	Straight	2 mm	41-70152
	Straight	3 mm	41-70153
	Straight	4 mm	41-70154

Titanium Straight Single-unit Components

	DESCRIPTION	REF #
	Closed Tray Single-unit Impression Coping	56-70021
	Open Tray Single-unit Impression Coping	56-70023
 	Single-unit Analog	56-70022
	Single-unit Digital Analog	30-70507
	Single-unit Healing Cap	55-70012
 	Titanium Single-unit Cylinder Temporary	41-70056
	Titanium Single-unit Cylinder Non-Engage	41-70016

*For use with the Closed Tray Multi-unit Impression Coping (56-72003).

Titanium Straight Single-unit Components cont.

	DESCRIPTION	REF #
	Fully Castable Single-unit Abutment	41-70057
	Fully Castable Single-unit Abutment Non-Engage	41-70017
	Titanium Single-unit Interface Coping	41-70059
	Single-unit Scan Abutment	30-70104

Straight Ball Abutments

	DESCRIPTION	CUFF	REF #
	Straight Ball Abutment	1 mm	40-70033
	Straight Ball Abutment	2 mm	40-70034
	Straight Ball Abutment	3 mm	40-70035

Straight Ball Components

	DESCRIPTION	REF #
	Ball Abutment Analog	56-70013
	Nylon Ball Abutment Retentive Cap Very Elastic - 350gr, Green	40-70131
	Nylon Ball Abutment Retentive Cap Extra Soft - 500gr, Yellow	40-70132
	Nylon Ball Abutment Retentive Cap Soft - 900gr, Pink	40-70023
	Nylon Ball Abutment Retentive Cap Standard - 1300gr, Transparent	40-70065
	Stainless Steel Ball Abutment Housing	40-70022

Paltop Equator Abutments

	DESCRIPTION	CUFF	REF #
	Straight	0.5 mm	40-70244
	Straight	1 mm	40-70245
	Straight	2 mm	40-70246
	Straight	3 mm	40-70247
	Straight	4 mm	40-70248
	Straight	5 mm	40-70249
	Straight	6 mm	40-70250

Paltop Equator Components

	DESCRIPTION	REF #
	Stainless Steel Housing	40-70287
	Retentive Cap Rigid- 2700gr, Violet	40-70282
	Retentive Cap Standard- 1800gr, Clear	40-70283
	Retentive Cap Soft- 1200 gr, Pink	40-70284
	Retentive Cap Extra Soft- 600gr, Yellow	40-70285
	Retentive Cap for Lab, Black	40-70286
	Smart Box Housing With Black Cap	40-70305
	Protective Disc, Clear	40-70288
	Protective Disc, Pink	40-70289
	Impression Coping (with Interchangeable Cap)	56-70042
	Analog	56-70041
	Core Tool	60-70161
	Retentive Cap Assortment Kit	40-70306
	Retentive Cap Smart Box Kit	40-70307

Angulation Corrective System (ACS)

	DESCRIPTION	CUFF	REF #
	Titanium Base Abutment Engage	1 mm	40-70231
	Titanium Base Abutment Engage	3 mm	40-70233
	Titanium Base Abutment Non-Engage	1 mm	40-70232
	Titanium Base Abutment Non-Engage	3 mm	40-70234

Angulation Corrective System (ACS) Components

	DESCRIPTION	REF #
	Screw NP	80-70038
	Manual Hex Driver for Angulation Corrective System	60-70152
	Hex Key for Angulation Corrective System	60-70153

NT Trading Components

	DESCRIPTION	REF #
	Scan Abutment from NT-Trading	30-70108
	Milling Blank Abutment from NT-Trading, Ø 11.95 mm	30-70235
	Digital Analog from NT-Trading	30-70510
	Mill Lab Analog for NT-Trading	30-70239

Screws


	DESCRIPTION	REF #
	Open Tray Impression Coping Screw NP, Slim	80-70033
	Open Tray Impression Coping Screw NP	80-70008
	Closed Tray Impression Coping Screw NP	80-70009
	Abutment Screw, M1.6 Direct NP	80-70010
	Angulation Corrective System Screw NP	80-70038

IH INTERNAL HEX CONNECTION

3.75 mm

4.2 mm

5 mm




STANDARD

PLATFORM

STANDARD

WORKFLOW



Advanced Implant

ADVANCED CLASSIC	DESCRIPTION	REF #
	Ø 3.75 X 8 mm	20-70017
	Ø 3.75 X 10 mm	20-70001
	Ø 3.75 X 11.5 mm	20-70002
	Ø 3.75 X 13 mm	20-70003
	Ø 3.75 X 16 mm	20-70004
	Ø 4.2 X 6 mm	20-70005
	Ø 4.2 X 8 mm	20-70006
	Ø 4.2 X 10 mm	20-70007
	Ø 4.2 X 11.5 mm	20-70008
	Ø 4.2 X 13 mm	20-70009
	Ø 4.2 X 16 mm	20-70010
	Ø 5 X 6 mm	20-70011
	Ø 5 X 8 mm	20-70012
	Ø 5 X 10 mm	20-70013
	Ø 5 X 11.5 mm	20-70014
	Ø 5 X 13 mm	20-70015
	Ø 5 X 16 mm	20-70016

Advanced + Implant

ADVANCED+ SLIGHTLY ACTIVE APEX	DESCRIPTION	REF #
	Ø 3.75 X 8 mm	20-70017P
	Ø 3.75 X 10 mm	20-70001P
	Ø 3.75 X 11.5 mm	20-70002P
	Ø 3.75 X 13 mm	20-70003P
	Ø 3.75 X 16 mm	20-70004P
	Ø 4.2 X 6 mm	20-70005P
	Ø 4.2 X 8 mm	20-70006P
	Ø 4.2 X 10 mm	20-70007P
	Ø 4.2 X 11.5 mm	20-70008P
	Ø 4.2 X 13 mm	20-70009P
	Ø 4.2 X 16 mm	20-70010P
	Ø 5 X 6 mm	20-70011P
	Ø 5 X 8 mm	20-70012P
	Ø 5 X 10 mm	20-70013P
	Ø 5 X 11.5 mm	20-70014P
	Ø 5 X 13 mm	20-70015P
	Ø 5 X 16 mm	20-70016P

Dynamic Implant

	DESCRIPTION	REF #
	Ø 3.75 X 8 mm	21-70017
	Ø 3.75 X 10 mm	21-70001
	Ø 3.75 X 11.5 mm	21-70002
	Ø 3.75 X 13 mm	21-70003
	Ø 3.75 X 16 mm	21-70004
	Ø 4.2 X 6 mm	21-70005
	Ø 4.2 X 8 mm	21-70006
	Ø 4.2 X 10 mm	21-70007
	Ø 4.2 X 11.5 mm	21-70008
	Ø 4.2 X 13 mm	21-70009
	Ø 4.2 X 16 mm	21-70010
	Ø 5 X 6 mm	21-70011
	Ø 5 X 8 mm	21-70012
	Ø 5 X 10 mm	21-70013
	Ø 5 X 11.5 mm	21-70014
	Ø 5 X 13 mm	21-70015
	Ø 5 X 16 mm	21-70016

DIVA by Paltop Dental Implant Set

	DESCRIPTION	REF #
	Ø 3.75 X 11.5 mm	23-70001S*
	Ø 3.75 X 13 mm	23-70002S*
	Ø 4.2 X 11.5 mm	23-70003S*
	Ø 4.2 X 13 mm	23-70004S*
	Ø 5 X 11.5 mm	23-70005S*
	Ø 5 X 13 mm	23-70006S*

*The DIVA by Paltop Dental Implant Set contains:
 DIVA implant with cover screw | Internal screw | Hex Key for DIVA Ø1mm, Manual | Syringe 10ml | Venflon | Cerasorb Paste.

PAI Implant

	DESCRIPTION	REF #
	Ø 3.75 X 8 mm	29-70017
	Ø 3.75 X 10 mm	29-70001
	Ø 3.75 X 11.5 mm	29-70002
	Ø 3.75 X 13 mm	29-70003
	Ø 3.75 X 16 mm	29-70004
	Ø 4.2 X 6 mm	29-70005
	Ø 4.2 X 8 mm	29-70006
	Ø 4.2 X 10 mm	29-70007
	Ø 4.2 X 11.5 mm	29-70008
	Ø 4.2 X 13 mm	29-70009
	Ø 4.2 X 16 mm	29-70010
	Ø 5 X 6 mm	29-70011
	Ø 5 X 8 mm	29-70012
	Ø 5 X 10 mm	29-70013
	Ø 5 X 11.5 mm	29-70014
	Ø 5 X 13 mm	29-70015
	Ø 5 X 16 mm	29-70016

PAI Total Coverage (TC) Implant

	DESCRIPTION	REF #
	Ø 3.75 X 8 mm	29-70017TC
	Ø 3.75 X 10 mm	29-70001TC
	Ø 3.75 X 11.5 mm	29-70002TC
	Ø 3.75 X 13 mm	29-70003TC
	Ø 3.75 X 16 mm	29-70004TC
	Ø 4.2 X 6 mm	29-70005TC
	Ø 4.2 X 8 mm	29-70006TC
	Ø 4.2 X 10 mm	29-70007TC
	Ø 4.2 X 11.5 mm	29-70008TC
	Ø 4.2 X 13 mm	29-70009TC
	Ø 4.2 X 16 mm	29-70010TC
	Ø 5 X 6 mm	29-70011TC
	Ø 5 X 8 mm	29-70012TC
	Ø 5 X 10 mm	29-70013TC
	Ø 5 X 11.5 mm	29-70014TC
	Ø 5 X 13 mm	29-70015TC
	Ø 5 X 16 mm	29-70016TC

Cover Screw

	DESCRIPTION	REF #
	SP	80-70001

Concave Healing Caps

	DESCRIPTION	CUFF	REF #
	Ø 4.5 mm	1 mm	55-70003
	Ø 4.5 mm	2 mm	55-70001
	Ø 4.5 mm	3 mm	55-70004
	Ø 4.5 mm	4 mm	55-70002
	Ø 4.5 mm	5 mm	55-70013
	Ø 4.5 mm	6 mm	55-70020
	Ø 4.5 mm	7 mm	55-70021
	Ø 5.5 mm	1 mm	55-70005
	Ø 5.5 mm	2 mm	55-70006
	Ø 5.5 mm	3 mm	55-70007
	Ø 5.5 mm	4 mm	55-70008
	Ø 5.5 mm	5 mm	55-70016

Straight Healing Caps

	DESCRIPTION	CUFF	REF #
	Ø 4.5 mm	2 mm	55-70025
	Ø 4.5 mm	3 mm	55-70026
	Ø 4.5 mm	5 mm	55-70027
	Ø 5.5 mm	2 mm	55-70034
	Ø 5.5 mm	3 mm	55-70035
	Ø 5.5 mm	5 mm	55-70036

Impression Copings

	DESCRIPTION	REF #
	Open Tray Impression Coping Ø 4.5 mm	56-70005
 	Open Tray Impression Coping, Slim	56-70032
	Open Tray Impression Coping, Short	56-70031
 	Closed Tray Impression Coping 15°, Angulated	56-70043
	Closed Tray Impression Coping 25°, Angulated	56-70044
	Closed Tray Impression Coping	56-70006
	Plastic Snap-On Closed Tray Impression Coping With Abutment	56-70034

Scan Abutment

	DESCRIPTION	REF #
	SP	30-70101

Implant Analogs

	DESCRIPTION	REF #
 	Implant Analog	56-70002
	Digital Analog	30-70501

Titanium Temporary Abutments

	DESCRIPTION	CUFF	REF #
	Engage	0.5 mm	45-70029
	Engage	2 mm	45-70051
	Engage	3 mm	45-70053
	Non-Engage	0.5 mm	45-70030
	Non-Engage	2 mm	45-70052
	Non-Engage	3 mm	45-70054

Titanium Temporary Abutments, Immediate

	DESCRIPTION	CUFF	REF #
	SP	1.5 mm	45-70033
	SP	3 mm	45-70035
	SP, Wide	3 mm	45-70050
	Hex Driver for Ratchet for Immediate Temporary Abutment		60-70244


PEEK Concave Temporary Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	45-70001
	Straight	2 mm	45-70002
	Straight	3 mm	45-70003
	Angulated 15°	1 mm	45-70006
	Angulated 15°	2 mm	45-70007
	Angulated 15°	3 mm	45-70008

Fully Castable Straight Abutment

	DESCRIPTION	REF #
	Engage	46-70013
	Non-Engage	46-70018

Castable Abutments

	DESCRIPTION	REF #
	Gold Base Engage	46-70020
	Gold Base Non-Engage	46-70022
	Cobalt Chrome Base Engage	40-70077
	Cobalt Chrome Base Non-Engage	40-70177

Titanium Straight Abutments

	DESCRIPTION	CUFF	REF #
	Ø 4.5 mm	1 mm	40-70001
	Ø 4.5 mm	2 mm	40-70002
	Ø 4.5 mm	3 mm	40-70003
	Ø 4.5 mm	4 mm	40-70012
	Anatomic Ø 4.5 mm	1 mm	40-70041
	Anatomic Ø 4.5 mm	2 mm	40-70042
	Anatomic Ø 4.5 mm	3 mm	40-70043
	Ø 5.5 mm	1 mm	40-70014
	Ø 5.5 mm	2 mm	40-70015
	Ø 5.5 mm	3 mm	40-70016
	Anatomic Concave Ø 4.5 mm	1 mm	40-70044
	Anatomic Concave Ø 4.5 mm	2 mm	40-70045
	Anatomic Concave Ø 4.5 mm	3 mm	40-70046
	Knife Edge Ø 3.8 mm		40-70013
	Slim Knife Edge Ø 3.8 mm		40-70115

Titanium Angulated Abutments

	DESCRIPTION	CUFF	REF #
	Anatomic 15°	1 mm	40-70006
	Anatomic 15°	2 mm	40-70007
	Anatomic 15°	3 mm	40-70008
	Anatomic Concave 15°	1 mm	40-70047
	Anatomic Concave 15°	2 mm	40-70048
	Anatomic Concave 15°	3 mm	40-70049
	Anatomic Concave 15°	4 mm	40-70074
	Anatomic 25°	1 mm	40-70009
	Anatomic 25°	2 mm	40-70010
	Anatomic 25°	3 mm	40-70011
	Anatomic Concave 25°	1 mm	40-70050
	Anatomic Concave 25°	2 mm	40-70051
	Anatomic Concave 25°	3 mm	40-70052
	Knife Edge 15°		40-70004
	Slim Knife Edge 15°		40-70116
	Knife Edge 25°		40-70005

Snap-On Abutment System (SAS)

	DESCRIPTION	LENGTH	CUFF	REF #
	Ø 4.5 mm	4 mm	1 mm	44-70001
	Ø 4.5 mm	4 mm	2 mm	44-70002
	Ø 4.5 mm	4 mm	3 mm	44-70003
	Ø 4.5 mm	4 mm	4 mm	44-70004
	Ø 4.5 mm	4 mm	5 mm	44-70005
	Ø 4.5 mm	5.5 mm	1 mm	44-70006
	Ø 4.5 mm	5.5 mm	2 mm	44-70007
	Ø 4.5 mm	5.5 mm	3 mm	44-70008
	Ø 4.5 mm	5.5 mm	4 mm	44-70009
	Ø 4.5 mm	5.5 mm	5 mm	44-70010
	Ø 6 mm	4 mm	1 mm	44-70011
	Ø 6 mm	4 mm	2 mm	44-70012
	Ø 6 mm	4 mm	3 mm	44-70013
	Ø 6 mm	4 mm	4 mm	44-70014
	Ø 6 mm	4 mm	5 mm	44-70015

Snap-On Abutment System (SAS) cont.

	DESCRIPTION	LENGTH	CUFF	REF #
	Ø 6 mm	5.5 mm	1 mm	44-70016
	Ø 6 mm	5.5 mm	2 mm	44-70017
	Ø 6 mm	5.5 mm	3 mm	44-70018
	Ø 6 mm	5.5 mm	4 mm	44-70019
	Ø 6 mm	5.5 mm	5 mm	44-70020

Titanium Base Abutments Round

	DESCRIPTION	LENGTH	CUFF	REF #
	Engage			40-70207
	Engage	4 mm	0.5 mm	40-70318
	Engage	4 mm	1 mm	40-70291
	Engage	4 mm	2 mm	40-70209
	Engage	4 mm	3 mm	
	Engage			40-70211
	Engage	6 mm	0.5 mm	40-70320
	Engage	6 mm	1 mm	40-70293
	Engage	6 mm	2 mm	40-70213
	Engage	6 mm	3 mm	
	Non-Engage			40-70208
	Non-Engage	4 mm	0.5 mm	40-70319
	Non-Engage	4 mm	1 mm	40-70292
	Non-Engage	4 mm	2 mm	40-70210
	Non-Engage	4 mm	3 mm	
	Non-Engage			40-70212
	Non-Engage	6 mm	0.5 mm	40-70321
	Non-Engage	6 mm	1 mm	40-70294
	Non-Engage	6 mm	2 mm	40-70214
	Non-Engage	6 mm	3 mm	

Titanium Base Abutments Round Anterior

	DESCRIPTION	LENGTH	CUFF	REF #
	Engage	6 mm	0.5 mm	40-70299
	Engage	6 mm	2 mm	40-70300
	Engage	6 mm	3 mm	40-70301

Titanium Multi-unit Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	42-70001
	Straight	2 mm	42-70002
	Straight	3 mm	42-70003
	Straight	4 mm	42-70004
	Straight	5 mm	42-70005
	Angulated 17°	2.5 mm	42-70006
	Angulated 17°	3.25 mm	42-70010
	Angulated 17°	4 mm	42-70007
	Angulated 30°	3.5 mm	42-70008
	Angulated 30°	4.5 mm	42-70009

Titanium Multi-unit Components

	DESCRIPTION	REF #
	Closed Tray Multi-unit Impression Copping	56-72003
	Open Tray Multi-unit Impression Copping	56-72002
	Plastic Snap-On Closed Tray Impression Copping	56-70049
	Plastic Snap-On Closed Tray Impression Copping, Ø4.5mm, (5 pack)	56-70062
	Multi-unit Analog	56-72001
	Multi-unit Digital Analog	30-70504
	Multi-unit Healing Cap, Tapered	55-72011
	Multi-unit Healing Cap	55-72010
	Titanium Multi-unit Cylinder H 10 mm	42-70020
	Titanium Multi-unit Cylinder H 15 mm	42-70030
	Fully Castable Multi-unit Abutment	42-70021
	Titanium Retentive Multi-unit Interface Copping	42-70028
	Multi-unit Scan Abutment	30-70103

Titanium Multi-unit Screws

	DESCRIPTION	REF #
	Multi-unit Screw UNF 1-72	80-70005
	Open Tray Multi-unit Impression Coping Screw UNF 1-72	80-70028
	Abutment Screw UNF 1-72 Direct	80-70002*

Titanium Straight Single-unit Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	41-70051
	Straight	2 mm	41-70052
	Straight	3 mm	41-70053
	Straight	4 mm	41-70054

Titanium Straight Single-unit Components

	DESCRIPTION	REF #
	Closed Tray Single-unit Impression Coping	56-70021
	Open Tray Single-unit Impression Coping	56-70023
 	Single-unit Analog	56-70022
	Single-unit Digital Analog	30-70507
	Single-unit Healing Cap	55-70012
 	Titanium Single-unit Cylinder Temporary	41-70056
	Titanium Single-unit Cylinder Non-Engage	41-70016
 	Fully Castable Single-unit Abutment	41-70057
	Fully Castable Single-unit Abutment Non-Engage	41-70017
	Titanium Single-unit Interface Coping	41-70059
	Single-unit Scan Abutment	30-70104

*For use with the Closed Tray Multi-unit Impression Coping (56-72003).

Straight Ball Abutments

	DESCRIPTION	CUFF	REF #
	Straight Ball Abutment	1 mm	40-70017
	Straight Ball Abutment	2 mm	40-70018
	Straight Ball Abutment	3 mm	40-70019
	Straight Ball Abutment	4 mm	40-70290

Straight Ball Components

	DESCRIPTION	REF #
	Ball Abutment Analog	56-70013
	Nylon Ball Abutment Retentive Cap Very Elastic - 350gr, Green	40-70131
	Nylon Ball Abutment Retentive Cap Extra Soft - 500gr, Yellow	40-70132
	Nylon Ball Abutment Retentive Cap Soft - 900gr, Pink	40-70023
	Nylon Ball Abutment Retentive Cap Standard - 1300gr, Transparent	40-70065
	Stainless Steel Ball Abutment Housing	40-70022

Paltop Equator Abutments

	DESCRIPTION	CUFF	REF #
	Straight	0.5 mm	40-70237
	Straight	1 mm	40-70238
	Straight	2 mm	40-70239
	Straight	3 mm	40-70240
	Straight	4 mm	40-70241
	Straight	5 mm	40-70242
	Straight	6 mm	40-70243

Paltop Equator Components

	DESCRIPTION	REF #
	Stainless Steel Housing	40-70287
	Retentive Cap Rigid- 2700gr, Violet	40-70282
	Retentive Cap Standard- 1800gr, Clear	40-70283
	Retentive Cap Soft- 1200 gr, Pink	40-70284

Paltop Equator Components cont.

	DESCRIPTION	REF #
	Retentive Cap Extra Soft- 600gr, Yellow	40-70285
	Retentive Cap for Lab, Black	40-70286
	Smart Box Housing With Black Cap	40-70305
	Protective Disc, Clear	40-70288
	Protective Disc, Pink	40-70289
	Impression Coping (with Interchangeable Cap)	56-70042
	Analog	56-70041
	Core Tool	60-70161
	Retentive Cap Assortment Kit	40-70306
	Retentive Cap Smart Box Kit	40-70307

Angulation Corrective System (ACS)

	DESCRIPTION	CUFF	REF #
	Titanium Base Abutment Engage	1 mm	40-70227
	Titanium Base Abutment Engage	3 mm	40-70229
	Titanium Base Abutment Non-Engage	1 mm	40-70228
	Titanium Base Abutment Non-Engage	3 mm	40-70230

Angulation Corrective System (ACS) Components

	DESCRIPTION	REF #
	Screw SP	80-70037
	Manual Hex Driver for Angulation Corrective System	60-70152
	Hex Key for Angulation Corrective System	60-70153


NT Trading Components


	DESCRIPTION	REF #
	Scan Abutment from NT-Trading	30-70107
	Milling Blank Abutment from NT-Trading, Ø 11.95 mm	30-70234
	Digital Analog from NT-Trading	30-70509
	Mill Lab Analog for NT-Trading	30-70238

Screws


	DESCRIPTION	REF #
	Open Tray Impression Coping Screw SP, Short	80-70031
	Open Tray Impression Coping Screw SP, Slim	80-70032
	Open Tray Impression Coping Screw UNF 1-72 SP	80-70003
	Closed Tray Impression Coping Screw UNF 1-72 SP	80-70004
	Abutment Screw UNF 1-72 Direct	80-70002
	Angulation Corrective System Screw SP	80-70037

IH INTERNAL HEX CONNECTION





WORKFLOW



Advanced + Implant

ADVANCED+ SLIGHTLY ACTIVE APEX	DESCRIPTION	REF #
	Ø 6 X 6 mm	20-70022P
	Ø 6 X 8 mm	20-70023P
	Ø 6 X 10 mm	20-70024P
	Ø 6 X 11.5 mm	20-70025P
	Ø 6 X 13 mm	20-70026P
	Ø 6 X 16 mm	20-70027P

Dynamic Implant

Dynamic	DESCRIPTION	REF #
	Ø 6 X 6 mm	21-70022
	Ø 6 X 8 mm	21-70023
	Ø 6 X 10 mm	21-70024
	Ø 6 X 11.5 mm	21-70025
	Ø 6 X 13 mm	21-70026
	Ø 6 X 16 mm	21-70027

PAI Implant

PAI	DESCRIPTION	REF #
	Ø 6 X 6 mm	29-70022
	Ø 6 X 8 mm	29-70023
	Ø 6 X 10 mm	29-70024
	Ø 6 X 11.5 mm	29-70025
	Ø 6 X 13 mm	29-70026
	Ø 6 X 16 mm	29-70027

PAI Total Coverage (TC) Implant

PAI	DESCRIPTION	REF #
	Ø 6 X 6 mm	29-70022TC
	Ø 6 X 8 mm	29-70023TC
	Ø 6 X 10 mm	29-70024TC
	Ø 6 X 11.5 mm	29-70025TC
	Ø 6 X 13 mm	29-70026TC
	Ø 6 X 16 mm	29-70027TC

Cover Screw

	DESCRIPTION	REF #
	WP	80-70029

Concave Healing Caps

	DESCRIPTION	CUFF	REF #
	Ø 6.5 mm	2 mm	55-70022
	Ø 6.5 mm	4 mm	55-70023
	Ø 6.5 mm	6 mm	55-70024

Straight Healing Caps

	DESCRIPTION	CUFF	REF #
	Ø 6.5 mm	2 mm	55-70031
	Ø 6.5 mm	4 mm	55-70032
	Ø 6.5 mm	6 mm	55-70033

Impression Copings

	DESCRIPTION	REF #
	Open Tray Impression Coping	56-70029
	Closed Tray Impression Coping	56-70030
	Plastic Snap-On Closed Tray Impression Coping With Abutment	56-70036


Scan Abutment

	DESCRIPTION	REF #
	WP	30-70106

Implant Analogs

	DESCRIPTION	REF #
	Implant Analog	56-70028
	Digital Analog	30-70506

Titanium Temporary Abutments

	DESCRIPTION	CUFF	REF #
	Engage	0.5 mm	45-70042
	Engage	2 mm	45-70059
	Engage	3 mm	45-70061
	Non-Engage	0.5 mm	45-70043
	Non-Engage	2 mm	45-70060
	Non-Engage	3 mm	45-70062

Titanium Temporary Abutments, Immediate

	DESCRIPTION	CUFF	REF #
	WP	1.5 mm	45-70048
	WP	3 mm	45-70049
	Hex Driver for Ratchet for Immediate Temporary Abutment		60-70244

Fully Castable Straight Abutment

	DESCRIPTION	REF #
	Engage	46-70033
	Non-Engage	46-70034

Castable Abutments

	DESCRIPTION	REF #
	Gold Base Engage	46-70041
	Gold Base Non-Engage	46-70042
	Cobalt Chrome Base Engage	40-70174
	Cobalt Chrome Base Non-Engage	40-70175

Titanium Straight Abutments

	DESCRIPTION	CUFF	REF #
	Anatomic Ø 6 mm	1 mm	40-70124
	Anatomic Ø 6 mm	2 mm	40-70125
	Anatomic Ø 6 mm	3 mm	40-70126
	Anatomic Ø 6 mm	4 mm	40-70127
	Anatomic Concave Ø 6 mm	1 mm	40-70120
	Anatomic Concave Ø 6 mm	2 mm	40-70121
	Anatomic Concave Ø 6 mm	3 mm	40-70122
	Anatomic Concave Ø 6 mm	4 mm	40-70123

Snap-On Abutment System (SAS)

	DESCRIPTION	LENGTH	CUFF	REF #
	Ø 6.0 mm	4.0 mm	1 mm	44-70041
	Ø 6.0 mm	4.0 mm	2 mm	44-70042
	Ø 6.0 mm	4.0 mm	3 mm	44-70043
	Ø 6.0 mm	5.5 mm	1 mm	44-70044
	Ø 6.0 mm	5.5 mm	2 mm	44-70045
	Ø 6.0 mm	5.5 mm	3 mm	44-70046
				


Titanium Base Abutments Round

	DESCRIPTION	LENGTH	CUFF	REF #
	Engage			40-70223
	Engage	4 mm	0.5 mm	40-70326
		4 mm	2 mm	
	Engage			40-70225
	Engage	6 mm	0.5 mm	40-70328
		6 mm	2 mm	
	Non-Engage			40-70224
	Non-Engage	4 mm	0.5 mm	40-70327
		4 mm	2 mm	
	Non-Engage			40-70226
	Non-Engage	6 mm	0.5 mm	40-70329
		6 mm	2 mm	

Titanium Multi-unit Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	42-70023
	Straight	2 mm	42-70024
	Straight	3 mm	42-70025
	Straight	4 mm	42-70026
	Angulated 17°	3 mm	42-70027

Titanium Multi-unit Components

	DESCRIPTION	REF #
	Closed Tray Multi-unit Impression Coping	56-72003
	Open Tray Multi-unit Impression Coping	56-72002
	Plastic Snap-On Closed Tray Impression Coping	56-70049
	Plastic Snap-On Closed Tray Impression Coping, Ø4.5mm, (5 pack)	56-70062
	Multi-unit Analog	56-72001
	Multi-unit Digital Analog	30-70504
	Multi-unit Healing Cap, Tapered	55-72011
	Multi-unit Healing Cap	55-72010
	Titanium Multi-unit Cylinder H 10 mm	42-70020
	Titanium Multi-unit Cylinder H 15 mm	42-70030
	Fully Castable Multi-unit Abutment	42-70021
	Titanium Retentive Multi-unit Interface Coping	42-70028
	Multi-unit Scan Abutment	30-70103

Titanium Multi-unit Screws

	DESCRIPTION	REF #
	Multi-unit Screw UNF 1-72	80-70005
	Open Tray Multi-unit Impression Coping Screw UNF 1-72	80-70028
	Abutment Screw UNF 1-72 Direct	80-70002*

Titanium Straight Single-unit Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	41-70159
	Straight	2 mm	41-70160
	Straight	3 mm	41-70161
	Straight	4 mm	41-70162

Titanium Straight Single-unit Components

	DESCRIPTION	REF #
	Closed Tray Single-unit Impression Coping	56-70021
	Open Tray Single-unit Impression Coping	56-70023
 	Single-unit Analog Single-unit Digital Analog	56-70022 30-70507
	Single-unit Healing Cap	55-70012
 	Titanium Single-unit Cylinder Temporary Titanium Single-unit Cylinder Non-Engage	41-70056 41-70016
 	Fully Castable Single-unit Abutment Fully Castable Single-unit Abutment Non-Engage	41-70057 41-70017
	Titanium Single-unit Interface Coping	41-70059
	Single-unit Scan Abutment	30-70104

*For use with the Closed Tray Multi-unit Impression Coping (56-72003).

Paltop Equator Abutments

	DESCRIPTION	CUFF	REF #
	Straight	0.5 mm	40-70251
	Straight	1 mm	40-70252
	Straight	2 mm	40-70253
	Straight	3 mm	40-70254
	Straight	4 mm	40-70255
	Straight	5 mm	40-70256
	Straight	6 mm	40-70257

Paltop Equator Components

	DESCRIPTION	REF #
	Stainless Steel Housing	40-70287
	Retentive Cap Rigid- 2700gr, Violet	40-70282
	Retentive Cap Standard- 1800gr, Clear	40-70283
	Retentive Cap Soft- 1200 gr, Pink	40-70284
	Retentive Cap Extra Soft- 600gr, Yellow	40-70285
	Retentive Cap for Lab, Black	40-70286
	Smart Box Housing With Black Cap	40-70305
	Protective Disc, Clear	40-70288
	Protective Disc, Pink	40-70289
	Impression Coping (with Interchangeable Cap)	56-70042
	Analog	56-70041
	Core Tool	60-70161
	Retentive Cap Assortment Kit	40-70306
	Retentive Cap Smart Box Kit	40-70307



NT Trading Components

	DESCRIPTION	REF #
	Scan Abutment from NT-Trading	30-70109
	Milling Blank Abutment from NT-Trading, Ø 11.95 mm	30-70236
	Digital Analog from NT-Trading	30-70511
	Mill Lab Analog for NT-Trading	30-70240

Screws

	DESCRIPTION	REF #
	Open Tray Impression Coping Screw WP	80-70036
	Closed Tray Impression Coping Screw UNF 1-72 WP	80-70030
	Abutment Screw UNF 1-72 Direct	80-70002

CC CONICAL CONNECTION




CONICAL
CONNECTION



CONICAL
ACTIVE

CONNECTION
WORKFLOW



PCA Implant



DESCRIPTION	REF #
Ø 3.25 X 10 mm	22-70018
Ø 3.25 X 11.5 mm	22-70019
Ø 3.25 X 13 mm	22-70020
Ø 3.25 X 16 mm	22-70021
Ø 3.75 X 8 mm	22-70017
Ø 3.75 X 10 mm	22-70001
Ø 3.75 X 11.5 mm	22-70002
Ø 3.75 X 13 mm	22-70003
Ø 3.75 X 16 mm	22-70004
Ø 4.2 X 6 mm	22-70005
Ø 4.2 X 8 mm	22-70006
Ø 4.2 X 10 mm	22-70007
Ø 4.2 X 11.5 mm	22-70008
Ø 4.2 X 13 mm	22-70009
Ø 4.2 X 16 mm	22-70010
Ø 5 X 6 mm	22-70011
Ø 5 X 8 mm	22-70012
Ø 5 X 10 mm	22-70013
Ø 5 X 11.5 mm	22-70014
Ø 5 X 13 mm	22-70015
Ø 5 X 16 mm	22-70016

Dynamic Conical Dental Implant



DESCRIPTION	REF #
Ø 3.25 X 10 mm	24-70018
Ø 3.25 X 11.5 mm	24-70019
Ø 3.25 X 13 mm	24-70020
Ø 3.25 X 16 mm	24-70021
Ø 3.75 X 8 mm	24-70017
Ø 3.75 X 10 mm	24-70001
Ø 3.75 X 11.5 mm	24-70002
Ø 3.75 X 13 mm	24-70003
Ø 3.75 X 16 mm	24-70004
Ø 4.2 X 6 mm	24-70005
Ø 4.2 X 8 mm	24-70006
Ø 4.2 X 10 mm	26-70007
Ø 4.2 X 11.5 mm	24-70008
Ø 4.2 X 13 mm	24-70009
Ø 4.2 X 16 mm	24-70010



Dynamic Conical Dental Implant cont.

	DESCRIPTION	REF #
	Ø 5 X 6 mm	24-70011
	Ø 5 X 8 mm	24-70012
	Ø 5 X 10 mm	24-70013
	Ø 5 X 11.5 mm	24-70014
	Ø 5 X 13 mm	24-70015
	Ø 5 X 16 mm	24-70016

Cover Screw

	DESCRIPTION	REF #
	Conical	80-70016

Concave Healing Caps

	DESCRIPTION	CUFF	REF #
	Ø 4.5 mm	2 mm	55-72042
	Ø 4.5 mm	3 mm	55-72043
	Ø 4.5 mm	4 mm	55-72044
	Ø 4.5 mm	5 mm	55-72045
	Ø 4.5 mm	7 mm	55-72047
	Ø 6 mm	2 mm	55-72062
	Ø 6 mm	3 mm	55-72063
	Ø 6 mm	5 mm	55-72065
	Ø 6 mm	7 mm	55-72067

Straight Healing Caps

 	DESCRIPTION	CUFF	REF #
	Ø 4.5 mm	2 mm	55-72002
	Ø 4.5 mm	3 mm	55-72003
	Ø 4.5 mm	5 mm	55-72005
	Ø 3 mm	2 mm	55-72006



Impression Copings

    	DESCRIPTION	REF #
	Open Tray Impression Coping Ø 4.5 mm	56-72045
	Open Tray Impression Coping Ø 4.5 mm, Slim	56-72046
	Closed Tray Impression Coping Ø 4.5 mm	56-72145
	Open Tray Impression Coping Ø 6 mm	56-72060
	Closed Tray Impression Coping Ø 6 mm	56-72160
	Plastic Snap-On Closed Tray Impression Coping With Abutment	56-72012

Scan Abutment

	DESCRIPTION	REF #
	Conical	30-70105

Implant Analogs

	DESCRIPTION	REF #
	Implant Analog	56-72011
	Digital Analog	30-70503

Titanium Temporary Abutments

 	DESCRIPTION	CUFF	REF #
	Engage		45-72031
	Engage	0.5 mm	45-72034
	Engage	2 mm	45-72036
		3 mm	
	Non-Engage		45-72032
	Non-Engage	0.5 mm	45-72035
	Non-Engage	2 mm	45-72037
		3 mm	

Titanium Temporary Abutments, Immediate

	DESCRIPTION	CUFF	REF #
	Conical	1.5 mm	45-72001
	Conical	3 mm	45-72003
	Conical, Wide	3 mm	45-72033
	Hex Driver for Ratchet for Immediate Temporary Abutment		60-70244

Castable Abutments

	DESCRIPTION	REF #
	Gold Base Engage	46-72001
	Gold Base Non-Engage	46-72002
	Cobalt Chrome Base Engage	40-72021
	Cobalt Chrome Base Non-Engage	40-72022

Titanium Straight Abutments

	DESCRIPTION	CUFF	REF #
	Anatomic Concave Ø 4.5 mm	1 mm	40-72041
	Anatomic Concave Ø 4.5 mm	2 mm	40-72042
	Anatomic Concave Ø 4.5 mm	3 mm	40-72043
	Anatomic Concave Ø 4.5 mm	4 mm	40-72044
	Anatomic Concave Ø 4.5 mm	5 mm	40-72045
	Anatomic Concave Ø 6 mm	1 mm	40-72061
	Anatomic Concave Ø 6 mm	2 mm	40-72062
	Anatomic Concave Ø 6 mm	3 mm	40-72063
	Anatomic Concave Ø 6 mm	4 mm	40-72064
	Anatomic Concave Ø 6 mm	5 mm	40-72065
	Knife Edge		40-72010

Titanium Angulated Abutments

	DESCRIPTION	CUFF	REF #
	Anatomic Concave 20° Ø4.5 mm	1 mm	40-72241
	Anatomic Concave 20° Ø4.5 mm	2 mm	40-72242
	Anatomic Concave 20° Ø4.5 mm	3 mm	40-72243
	Anatomic Concave 20° Ø4.5 mm	4 mm	40-72244
	Anatomic Concave 20° Ø4.5 mm	5 mm	40-72245
	Knife Edge 20°		40-72020

Snap-On Abutment System (SAS)

	DESCRIPTION	LENGTH	CUFF	REF #
	Ø 4.5 mm	4 mm	1 mm	44-72441
	Ø 4.5 mm	4 mm	2 mm	44-72442
	Ø 4.5 mm	4 mm	3 mm	44-72443
	Ø 4.5 mm	4 mm	4 mm	44-72444
	Ø 4.5 mm	4 mm	5 mm	44-72445
	Ø 6 mm	4 mm	1 mm	44-72461
	Ø 6 mm	4 mm	2 mm	44-72462
	Ø 6 mm	4 mm	3 mm	44-72463
	Ø 6 mm	4 mm	4 mm	44-72464
	Ø 6 mm	4 mm	5 mm	44-72465
	Ø 4.5 mm	5.5 mm	1 mm	44-72541
	Ø 4.5 mm	5.5 mm	2 mm	44-72542
	Ø 4.5 mm	5.5 mm	3 mm	44-72543
	Ø 4.5 mm	5.5 mm	4 mm	44-72544
	Ø 4.5 mm	5.5 mm	5 mm	44-72545
	Ø 6 mm	5.5 mm	1 mm	44-72561
	Ø 6 mm	5.5 mm	2 mm	44-72562
	Ø 6 mm	5.5 mm	3 mm	44-72563
	Ø 6 mm	5.5 mm	4 mm	44-72564
	Ø 6 mm	5.5 mm	5 mm	44-72565

Titanium Base Abutments Round

	DESCRIPTION	LENGTH	CUFF	REF #
	Engage	4 mm	0.5 mm	40-72046
	Engage	4 mm	1 mm	40-72088
	Engage	4 mm	2 mm	40-72084
	Engage	4 mm	3 mm	40-72048
	Engage	6 mm	0.5 mm	40-72050
	Engage	6 mm	1 mm	40-72090
	Engage	6 mm	2 mm	40-72086
	Engage	6 mm	3 mm	40-72052
	Non-Engage	4 mm	0.5 mm	40-72047
	Non-Engage	4 mm	1 mm	40-72089
	Non-Engage	4 mm	2 mm	40-72085
	Non-Engage	4 mm	3 mm	40-72049
	Non-Engage	6 mm	0.5 mm	40-72051
	Non-Engage	6 mm	1 mm	40-72091
	Non-Engage	6 mm	2 mm	40-72087
	Non-Engage	6 mm	3 mm	40-72053


Titanium Base Abutments Round Anterior

	DESCRIPTION	LENGTH	CUFF	REF #
	Engage	6 mm	0.5 mm	40-72081
	Engage	6 mm	2 mm	40-72082
	Engage	6 mm	3 mm	40-72083

Titanium Multi-unit Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	42-72041
	Straight	2 mm	42-72042
	Straight	3 mm	42-72043
	Straight	4 mm	42-72044
	Straight	5 mm	42-72045
	Angulated 17°	3 mm	42-72017
	Angulated 17°	4.5 mm	42-72018
	Angulated 30°	4 mm	42-72030

Titanium Multi-unit Components

	DESCRIPTION	REF #
	Closed Tray Multi-unit Impression Coping	56-72003
	Open Tray Multi-unit Impression Coping	56-72002
	Plastic Snap-On Closed Tray Impression Coping	56-70049
	Plastic Snap-On Closed Tray Impression Coping, Ø4.5mm, (5 pack)	56-70062
	Multi-unit Analog	56-72001
	Multi-unit Digital Analog	30-70504
	Multi-unit Healing Cap, Tapered	55-72011
	Multi-unit Healing Cap	55-72010
	Titanium Multi-unit Cylinder H 10 mm	42-70020
	Titanium Multi-unit Cylinder H 15 mm	42-70030
	Fully Castable Multi-unit Abutment	42-70021
	Titanium Retentive Multi-unit Interface Coping	42-70028
	Multi-unit Scan Abutment	30-70103

Titanium Multi-unit Screws

	DESCRIPTION	REF #
	Multi-unit Screw UNF 1-72	80-70005
	Open Tray Multi-unit Impression Coping Screw UNF 1-72	80-70028
	Abutment Screw UNF 1-72 Direct	80-70002*

*For use with the Closed Tray Multi-unit Impression Coping (56-72003)..

Titanium Straight Single-unit Abutments

	DESCRIPTION	CUFF	REF #
	Straight	1 mm	43-72041
	Straight	2 mm	43-72042
	Straight	3 mm	43-72043
	Straight	4 mm	43-72044
	Straight	5 mm	43-72045

Titanium Straight Single-unit Components

	DESCRIPTION	REF #
	Closed Tray Single-unit Impression Coping	56-70021
	Open Tray Single-unit Impression Coping	56-70023
 	Single-unit Analog Single-unit Digital Analog	56-70022 30-70507
	Single-unit Healing Cap	55-70012
 	Titanium Single-unit Cylinder Temporary Titanium Single-unit Cylinder Non-Engage	41-70056 41-70016
 	Fully Castable Single-unit Abutment Fully Castable Single-unit Abutment Non-Engage	41-70057 41-70017
	Titanium Single-unit Interface Coping	41-70059
	Single-unit Scan Abutment	30-70104

Straight Ball Abutments

	DESCRIPTION	CUFF	REF #
	Straight Ball Abutment	1 mm	40-72001
	Straight Ball Abutment	2 mm	40-72002
	Straight Ball Abutment	3 mm	40-72003
	Straight Ball Abutment	4 mm	40-72004
	Straight Ball Abutment	5 mm	40-72005

Straight Ball Components

	DESCRIPTION	REF #
	Ball Abutment Analog	56-70013
	Nylon Ball Abutment Retentive Cap Very Elastic - 350gr, Green	40-70131
	Nylon Ball Abutment Retentive Cap Extra Soft - 500gr, Yellow	40-70132
	Nylon Ball Abutment Retentive Cap Soft - 900gr, Pink	40-70023
	Nylon Ball Abutment Retentive Cap Standard - 1300gr, Transparent	40-70065
	Stainless Steel Ball Abutment Housing	40-70022


Paltop Equator Abutments

	DESCRIPTION	CUFF	REF #
	Straight	0.5 mm	40-72066
	Straight	1 mm	40-72067
	Straight	2 mm	40-72068
	Straight	3 mm	40-72069
	Straight	4 mm	40-72070
	Straight	5 mm	40-72071
	Straight	6 mm	40-72072

Paltop Equator Components

	DESCRIPTION	REF #
	Stainless Steel Housing	40-70287
	Retentive Cap Rigid- 2700gr, Violet	40-70282
	Retentive Cap Standard- 1800gr, Clear	40-70283
	Retentive Cap Soft- 1200 gr, Pink	40-70284
	Retentive Cap Extra Soft- 600gr, Yellow	40-70285
	Retentive Cap for Lab, Black	40-70286
	Smart Box Housing With Black Cap	40-70305
	Protective Disc, Clear	40-70288
	Protective Disc, Pink	40-70289

Paltop Equator Components cont.

	DESCRIPTION	REF #
	Impression Coping (with Interchangeable Cap)	56-70042
	Analog	56-70041
	Core Tool	60-70161
	Retentive Cap Assortment Kit	40-70306
	Retentive Cap Smart Box Kit	40-70307

Angulated Corrective System (ACS)

	DESCRIPTION	CUFF	REF #
	Titanium Base Abutment Engage	1 mm	40-72054
	Titanium Base Abutment Engage	3 mm	40-72056
	Titanium Base Abutment Non-Engage	1 mm	40-72055
	Titanium Base Abutment Non-Engage	3 mm	40-72057

Angulated Corrective System (ACS) Components

	DESCRIPTION	REF #
	Angulation Corrective System Screw Conical	80-72001
	Manual Hex Driver for Angulation Corrective System	60-70152
	Hex Key for Angulation Corrective System	60-70153

NT Trading Components

	DESCRIPTION	REF #
	Scan Abutment from NT-Trading	30-70110
	Milling Blank Abutment from NT-Trading, Ø 11.95 mm	30-70237
	Digital Analog from NT-Trading	30-70512
	Mill Lab Analog for NT-Trading	30-70241

Screws

	DESCRIPTION	REF #
	Open Tray Impression Coping Screw Conical	80-70012
	Closed Tray Impression Coping Screw Conical	80-70035
	Abutment Screw M1.6 Direct Conical	80-72002

IH INTERNAL HEX CONNECTION

CC CONICAL CONNECTION




LOCATOR[®]

PROSTHETIC
COMPONENTS

LOCATOR®

WORKFLOW



LOCATOR® Abutments


	DESCRIPTION	CUFF	REF #
	NP	1 mm	40-70088
	NP	2 mm	40-70089
	NP	3 mm	40-70090
	NP	4 mm	40-70091
	NP	5 mm	40-70092
	NP	6 mm	40-70093
	SP	1 mm	40-70082
	SP	2 mm	40-70083
	SP	3 mm	40-70084
	SP	4 mm	40-70085
	SP	5 mm	40-70086
	SP	6 mm	40-70087
	Conical	1 mm	40-72011
	Conical	2 mm	40-72012
	Conical	3 mm	40-72013
	Conical	4 mm	40-72014
	Conical	5 mm	40-72015
	Conical	6 mm	40-72016

LOCATOR® Components

	DESCRIPTION	REF #
	LOCATOR® Male Processing Package	40-70081
	LOCATOR® Male Processing Package Extended Range	40-70112
	LOCATOR® Impression Coping	56-70026
	LOCATOR® Female Analog, Ø 4 mm	56-70027
	Core Tool for LOCATOR®	60-70133

IH INTERNAL HEX CONNECTION

CC CONICAL CONNECTION




DIGITAL

PROSTHETIC
COMPONENTS


DIGITAL


ABUTMENT LEVEL
WORKFLOW




DIGITAL

IMPLANT LEVEL
WORKFLOW





CUSTOM MILLED SOLUTIONS
WORKFLOW



Scan Abutments

	DESCRIPTION	REF #
	Multi-unit Scan Abutment	30-70103
	Single-unit Scan Abutment	30-70104
	Scan Abutment NP	30-70102
	Scan Abutment SP	30-70101
	Scan Abutment WP	30-70106
	Scan Abutment Conical	30-70105

Digital Analogs

	DESCRIPTION	REF #
	Multi-unit Digital Analog	30-70504
	Single-unit Digital Analog	30-70507
	Digital Analog NP	30-70502
	Digital Analog SP	30-70501
	Digital Analog WP	30-70506
	Digital Analog Conical	30-70503


Interface Copings

	DESCRIPTION	REF #
	Titanium Retentive Multi-unit Interface Coping	42-70028
	Titanium Single-unit Interface Coping	41-70059

Titanium Base Abutments Round

	DESCRIPTION	LENGTH	CUFF	REF #
	Int. Hex NP Engage	4 mm	0.5 mm	40-70215
	Int. Hex NP Engage	4 mm	1 mm	40-70322
	Int. Hex NP Engage	4 mm	2 mm	40-70295
	Int. Hex NP Engage	4 mm	3 mm	40-70217
	Int. Hex NP Non-Engage	4 mm	0.5 mm	40-70216
	Int. Hex NP Non-Engage	4 mm	1 mm	40-70323
	Int. Hex NP Non-Engage	4 mm	2 mm	40-70296
	Int. Hex NP Non-Engage	4 mm	3 mm	40-70218
	Int. Hex NP Engage	6 mm	0.5 mm	40-70219
	Int. Hex NP Engage	6 mm	1 mm	40-70324
	Int. Hex NP Engage	6 mm	2 mm	40-70297
	Int. Hex NP Engage	6 mm	3 mm	40-70221
	Int. Hex NP Non-Engage	6 mm	0.5 mm	40-70220
	Int. Hex NP Non-Engage	6 mm	1 mm	40-70325
	Int. Hex NP Non-Engage	6 mm	2 mm	40-70298
	Int. Hex NP Non-Engage	6 mm	3 mm	40-70222




Titanium Base Abutments Round cont.

	DESCRIPTION	LENGTH	CUFF	REF #
	Int. Hex SP Engage	4 mm	0.5 mm	40-70207
	Int. Hex SP Engage	4 mm	1 mm	40-70318
	Int. Hex SP Engage	4 mm	2 mm	40-70291
	Int. Hex SP Engage	4 mm	3 mm	40-70209
	Int. Hex SP Non-Engage	4 mm	0.5 mm	40-70208
	Int. Hex SP Non-Engage	4 mm	1 mm	40-70319
	Int. Hex SP Non-Engage	4 mm	2 mm	40-70292
	Int. Hex SP Non-Engage	4 mm	3 mm	40-70210
	Int. Hex SP Engage	6 mm	0.5 mm	40-70211
	Int. Hex SP Engage	6 mm	1 mm	40-70320
	Int. Hex SP Engage	6 mm	2 mm	40-70293
	Int. Hex SP Engage	6 mm	3 mm	40-70213
	Int. Hex SP Non-Engage	6 mm	0.5 mm	40-70212
	Int. Hex SP Non-Engage	6 mm	1 mm	40-70321
	Int. Hex SP Non-Engage	6 mm	2 mm	40-70294
	Int. Hex SP Non-Engage	6 mm	3 mm	40-70214
	Int. Hex WP Engage	4 mm	0.5 mm	40-70223
	Int. Hex WP Engage	4 mm	2 mm	40-70326
	Int. Hex WP Non-Engage	4 mm	0.5 mm	40-70224
	Int. Hex WP Non-Engage	4 mm	2 mm	40-70327
	Int. Hex WP Engage	6 mm	0.5 mm	40-70225
	Int. Hex WP Engage	6 mm	2 mm	40-70328
	Int. Hex WP Non-Engage	6 mm	0.5 mm	40-70226
	Int. Hex WP Non-Engage	6 mm	2 mm	40-70329
	Conical Engage	4 mm	0.5 mm	40-72046
	Conical Engage	4 mm	1 mm	40-72088
	Conical Engage	4 mm	2 mm	40-72084
	Conical Engage	4 mm	3 mm	40-72048
	Conical Non-Engage	4 mm	0.5 mm	40-72047
	Conical Non-Engage	4 mm	1 mm	40-72089
	Conical Non-Engage	4 mm	2 mm	40-72085
	Conical Non-Engage	4 mm	3 mm	40-72049
	Conical Engage	6 mm	0.5 mm	40-72050
	Conical Engage	6 mm	1 mm	40-72090
	Conical Engage	6 mm	2 mm	40-72086
	Conical Engage	6 mm	3 mm	40-72052
	Conical Non-Engage	6 mm	0.5 mm	40-72051
	Conical Non-Engage	6 mm	1 mm	40-72091
	Conical Non-Engage	6 mm	2 mm	40-72087
	Conical Non-Engage	6 mm	3 mm	40-72053

Titanium Base Abutments Round Anterior

	DESCRIPTION	LENGTH	CUFF	REF #
	Int. Hex NP Engage	6 mm	0.5 mm	40-70302
	Int. Hex NP Engage	6 mm	2 mm	40-70303
	Int. Hex NP Engage	6 mm	3 mm	40-70304
	Int. Hex SP Engage	6 mm	0.5 mm	40-70299
	Int. Hex SP Engage	6 mm	2 mm	40-70300
	Int. Hex SP Engage	6 mm	3 mm	40-70301
	Conical Engage	6 mm	0.5 mm	40-72081
	Conical Engage	6 mm	2 mm	40-72082
	Conical Engage	6 mm	3 mm	40-72083

Titanium Temporary Abutments

	DESCRIPTION	CUFF	REF #
	Int. Hex NP Engage	0.5 mm	45-70031
	Int. Hex NP Engage	2 mm	45-70055
	Int. Hex NP Engage	3 mm	45-70057
	Int. Hex NP Non-Engage	0.5 mm	45-70032
	Int. Hex NP Non-Engage	2 mm	45-70056
	Int. Hex NP Non-Engage	3 mm	45-70058
	Int. Hex SP Engage	0.5 mm	45-70029
	Int. Hex SP Engage	2 mm	45-70051
	Int. Hex SP Engage	3 mm	45-70053
	Int. Hex SP Non-Engage	0.5 mm	45-70030
	Int. Hex SP Non-Engage	2 mm	45-70052
	Int. Hex SP Non-Engage	3 mm	45-70054
	Int. Hex WP Engage	0.5 mm	45-70042
	Int. Hex WP Engage	2 mm	45-70059
	Int. Hex WP Engage	3 mm	45-70061
	Int. Hex WP Non-Engage	0.5 mm	45-70043
	Int. Hex WP Non-Engage	2 mm	45-70060
	Int. Hex WP Non-Engage	3 mm	45-70062

Angulation Corrective System (ACS)

	DESCRIPTION	CUFF	REF #
	Titanium Base Abutment Int. Hex NP Engage	1 mm	40-70231
	Titanium Base Abutment Int. Hex NP Engage	3 mm	40-70233
	Titanium Base Abutment Int. Hex NP Non-Engage	1 mm	40-70232
	Titanium Base Abutment Int. Hex NP Non-Engage	3 mm	40-70234
	Titanium Base Abutment Int. Hex SP Engage	1 mm	40-70227
	Titanium Base Abutment Int. Hex SP Engage	3 mm	40-70229
	Titanium Base Abutment Int. Hex SP Non-Engage	1 mm	40-70228
	Titanium Base Abutment Int. Hex SP Non-Engage	3 mm	40-70230
	Titanium Base Abutment Conical Engage	1 mm	40-72054
	Titanium Base Abutment Conical Engage	3 mm	40-72056
	Titanium Base Abutment Conical Non-Engage	1 mm	40-72055
	Titanium Base Abutment Conical Non-Engage	3 mm	40-72057

Angulation Corrective System (ACS) Components

	DESCRIPTION	REF #
	Angulation Corrective System Screw NP	80-70038
	Angulation Corrective System Screw SP	80-70037
	Angulation Corrective System Screw Conical	80-72001
	Manual Hex Driver for Angulation Corrective System	60-70152
	Hex Key for Angulation Corrective System	60-70153

NT Trading Components

	DESCRIPTION	REF #
	Scan Abutment from NT-Trading Internal Hex. NP	30-70108
	Milling Blank Abutment from NT-Trading Internal Hex. NP Ø11.95mm	30-70235
	Digital Analog from NT-Trading NP	30-70510
	Mill Lab Analog for NT Trading NP	30-70239
	Scan Abutment from NT-Trading Internal Hex. SP	30-70107
	Milling Blank Abutment from NT-Trading Internal Hex. SP Ø11.95mm	30-70234
	Digital Analog from NT-Trading SP	30-70509
	Mill Lab Analog for NT Trading SP	30-70238

NT Trading cont.

	DESCRIPTION	REF #
	Scan Abutment from NT-Trading Internal Hex. WP	30-70109
	Milling Blank Abutment from NT-Trading Internal Hex. WP Ø 11.95mm	30-70236
	Digital Analog from NT-Trading WP	30-70511
	Mill Lab Analog for NT-Trading WP	30-70240
	Scan Abutment from NT-Trading Conical	30-70110
	Milling Blank Abutment from NT-Trading Conical Ø11.95mm	30-70237
	Digital Analog from NT-Trading Conical	30-70512
	Mill Lab Analog for NT-Trading Conical	30-70241



SURGICAL

KITS

PREMIUM

KIT





Premium Surgical Kit

DESCRIPTION	QUANTITY	REF #
Step Drill Ø 1.5/2 mm, Short		60-70047
Step Drill Ø 2/2.4 mm, Short		60-70049
Step Drill Ø 2.4/3.2 mm, Short		60-70051
Step Drill Ø 3.2/3.8 mm, Short		60-70053
Step Drill Ø 3.8/4.3 mm, Short		60-70120
Pilot Drill Spade		60-70020
Step Drill Ø 1.5/2 mm, Long		60-70048
Step Drill Ø 2/2.4 mm, Long		60-70050
Step Drill Ø 2.4/3.2 mm, Long		60-70052
Step Drill Ø 3.2/3.8 mm, Long		60-70054
Step Drill Ø 3.8/4.3 mm, Long		60-70121
Final Drill for Ø 3 mm Implant, Long		60-70147
Final Drill for Ø 3.25 mm Implant, Long		60-70035
Final Drill for Ø 3.75 mm Implant, Long		60-70037
Final Drill for Ø 4.2 mm Implant, Long		60-70039
Final Drill for Ø 5 mm Implant, Long		60-70041
Final Drill for Ø 6 mm Implant, Long		60-70098
Countersink for Ø 3.25 mm & Ø 3.75 mm Implants		60-70145
Countersink for Ø 4.2 mm Implant		60-70144
Countersink for Ø 5 mm & Ø 6 mm Implants		60-70143
Key Ratchet Adapter	2	60-70100
Drill Extender		60-70019

Premium Surgical Kit cont.

DESCRIPTION	QUANTITY	REF #
Hex Key Ø 1.25 mm L 15 mm, Long		60-70101
Hex Key Ø 2 mm, Long		60-70105
Hex Key Ø 2.5 mm, Long		60-70103
Hex Key Ø 1.25 mm L 9 mm, Short		60-70102
Hex Key Ø 2 mm, Short		60-70106
Hex Key Ø 2.5 mm, Short		60-70104
Parallel Pin Ø 1.5/2 mm	4	60-70015
Manual Hex Driver Ø 1.25 mm L 14 mm, Standard		60-70006
Torque Ratchet With Bending beam		60-70079
Drill Depth Probe		60-70073
Triangular Drill Stop Kit		60-70230
Hand Piece Straight		60-70014
Surgical Drilling Protocol		85-70041

Step Drills

	DESCRIPTION	REF #
	Step Drill Ø 1.5/2 mm, Short	60-70047
	Step Drill Ø 2/2.4 mm, Short	60-70049
	Step Drill Ø 2.4/3.2 mm, Short	60-70051
	Step Drill Ø 3.2/3.8 mm, Short	60-70053
	Step Drill Ø 3.8/4.3 mm, Short	60-70120
	Step Drill Ø 1.5/2 mm, Long	60-70048
	Step Drill Ø 2/2.4 mm, Long	60-70050
	Step Drill Ø 2.4/3.2 mm, Long	60-70052
	Step Drill Ø 3.2/3.8 mm, Long	60-70054
	Step Drill Ø 3.8/4.3 mm, Long	60-70121


Pilot Drill Spade

	DESCRIPTION	REF #
	Pilot Drill Spade	60-70020

Final Drills

	DESCRIPTION	REF #
	Final Drill for Ø 3 mm Implant, Long	60-70147
	Final Drill for Ø 3.25 mm Implant, Long	60-70035
	Final Drill for Ø 3.75 mm Implant, Long	60-70037
	Final Drill for Ø 4.2 mm Implant, Long	60-70039
	Final Drill for Ø 5 mm Implant, Long	60-70041
	Final Drill for Ø 6 mm Implant, Long	60-70098

Countersinks

	DESCRIPTION	REF #
	Countersink for Ø 3.25 mm & Ø 3.75 mm Implants	60-70145
	Countersink for Ø 4.2 mm Implant	60-70144
	Countersink for Ø 5 mm & Ø 6 mm Implants	60-70143

Key Ratchet Adapter

	DESCRIPTION	REF #
	Key Ratchet Adapter	60-70100

Drill Extender

	DESCRIPTION	REF #
	Drill Extender	60-70019


Hex Keys

	DESCRIPTION	REF #
	Hex Key Ø 1.25 mm L 15 mm, Long	60-70101
	Hex Key Ø 2 mm, Long	60-70105
	Hex Key Ø 2.5 mm, Long	60-70103
	Hex Key Ø 1.25 mm L 9 mm, Short	60-70102
	Hex Key Ø 2 mm, Short	60-70106
	Hex Key Ø 2.5 mm, Short	60-70104

Parallel Pins

	DESCRIPTION	REF #
	Parallel Pin Ø 1.5/2 mm	60-70015

Manual Hex Driver

	DESCRIPTION	REF #
	Manual Hex Driver Ø 1.25 mm L 14 mm, Standard	60-70006

Torque Ratchet

	DESCRIPTION	REF #
	Torque Ratchet With Bending beam	60-70079


Drill Depth Probe

	DESCRIPTION	REF #
	Drill Depth Probe	60-70073


Hand Piece

DESCRIPTION		REF #
	Hand Piece Straight	60-70014

Triangular Drill Stop Kit

DESCRIPTION		REF #
	Triangular Drill Stop Kit	60-70230*

**See pg. 99 for product details.*



KIT



Fully Guided Surgical Kit

30-70021

DESCRIPTION	REF #
Contra Angle NSK	30-50324
Digital Guidance Sleeve (DGS) for NSK, Narrow	30-70616
Digital Guidance Sleeve (DGS) for NSK	30-70600
Fully Guided Drill Ø 2 mm, Short	30-70305
Tissue Punch Ø 4 mm	60-70140
Tissue Punch Ø 5 mm	60-70141
Contra Angle Adapter	30-70618
Fully Guided Drill Ø 2 mm L 20 mm	30-70306
Fully Guided Drill Ø 3.25 mm L 20 mm	30-70307
Fully Guided Drill Ø 3.75 mm L 20 mm	30-70308
Fully Guided Drill Ø 4.2 mm L 20 mm	30-70309
Fully Guided Drill Ø 5 mm L 20 mm	30-70310
Fully Guided Drill Ø 2 mm L 25 mm	30-70312
Fully Guided Drill Ø 3.25 mm L 25 mm	30-70313
Fully Guided Drill Ø 3.75 mm L 25 mm	30-70314
Fully Guided Drill Ø 4.2 mm L 25 mm	30-70315
Fully Guided Drill Ø 5 mm L 25 mm	30-70316
Fully Guided Drill Ø 2 mm L 30 mm	30-70318
Fully Guided Drill Ø 3.25 mm L30 mm	30-70319
Fully Guided Drill Ø 3.75 mm L 30 mm	30-70320
Fully Guided Drill Ø 4.2 mm L 30 mm	30-70321
Fully Guided Drill Ø 5 mm L 30 mm	30-70322

Fully Guided Surgical Kit cont.

30-70021

DESCRIPTION	REF #
Implant Insertion Key NP + Conical, Narrow, Short	30-70640
Implant Insertion Key NP + Conical, Narrow, Medium	30-70641
Implant Insertion Key NP + Conical, Narrow, Long	30-70642
Implant Insertion Key SP, Narrow, Short	30-70622
Implant Insertion Key SP, Narrow, Medium	30-70632
Implant Insertion Key SP, Narrow	30-70602
Implant Insertion Key NP + Conical, Wide, Short	30-70643
Implant Insertion Key NP + Conical, Wide, Medium	30-70644
Implant Insertion Key NP + Conical, Wide, Long	30-70645
Implant Insertion Key SP, Wide, Short	30-70624
Implant Insertion Key SP, Wide, Medium	30-70634
Implant Insertion Key SP, Wide	30-70604
Manual Hex Driver Ø 1.25 mm L 14 mm, Standard	60-70006
Fully Guided Drill for Pin Fixation Ø 2 mm L 20 mm	30-70333
Fixation Pin L 20 mm (4 Pack)	30-70516
Torque Ratchet With Bending beam	60-70079
Fully Guided System Quick Start Brochure	85-70042

Contra Angle NSK

DESCRIPTION	REF #
 Contra Angle NSK	30-50324

Digital Guidance Sleeves (DGS)

DESCRIPTION	REF #
 Digital Guidance Sleeve (DGS) for NSK, Narrow	30-70616
Digital Guidance Sleeve (DGS) for NSK	30-70600

Short Pilot Drill

DESCRIPTION	REF #
 Fully Guided Drill Ø 2 mm, Short	30-70305

Tissue Punches

DESCRIPTION	REF #
 Tissue Punch Ø 4 mm	60-70140
Tissue Punch Ø 5 mm	60-70141

Contra Angle Adaptor

DESCRIPTION	REF #
 Contra Angle Adaptor	30-70618


Fully Guided Drills

DESCRIPTION	REF #
 Fully Guided Drill Ø 2 mm L 20 mm	30-70306
 Fully Guided Drill Ø 3.25 mm L 20 mm	30-70307
 Fully Guided Drill Ø 3.75 mm L 20 mm	30-70308
 Fully Guided Drill Ø 4.2 mm L 20 mm	30-70309
 Fully Guided Drill Ø 5 mm L 20 mm	30-70310
Fully Guided Drill Ø 2 mm L 25 mm	30-70312
Fully Guided Drill Ø 3.25 mm L 25 mm	30-70313
Fully Guided Drill Ø 3.75 mm L 25 mm	30-70314
Fully Guided Drill Ø 4.2 mm L 25 mm	30-70315
Fully Guided Drill Ø 5 mm L 25 mm	30-70316
Fully Guided Drill Ø 2 mm L 30 mm	30-70318
Fully Guided Drill Ø 3.25 mm L 30 mm	30-70319
Fully Guided Drill Ø 3.75 mm L 30 mm	30-70320
Fully Guided Drill Ø 4.2 mm L 30 mm	30-70321
Fully Guided Drill Ø 5 mm L 30 mm	30-70322

Implant Insertion Keys

DESCRIPTION	REF #
 NP + Conical, Narrow, Short	30-70640
 NP + Conical, Narrow, Medium	30-70641
 NP + Conical, Narrow, Long	30-70642
SP, Narrow, Short	30-70622
SP, Narrow, Medium	30-70632
SP, Narrow	30-70602
NP + Conical, Wide, Short	30-70643
NP + Conical, Wide, Medium	30-70644
NP + Conical, Wide, Long	30-70645
SP, Wide, Short	30-70624
SP, Wide, Medium	30-70634
SP, Wide	30-70604

Manual Hex Driver

DESCRIPTION	REF #
 Manual Hex Driver Ø 1.25 mm L 14 mm, Standard	60-70006

Fully Guided Drill for Pin Fixation

DESCRIPTION	REF #
 Fully Guided Drill for Pin Fixation Ø 2 mm L 20 mm	30-70333

Fixation Pin

DESCRIPTION	REF #
 Fixation Pin L 20 mm (4 Pack)	30-70516

Torque Ratchet

DESCRIPTION	REF #
 Torque Ratchet With Bending beam	60-70079

Guided Countersink*

	DESCRIPTION	REF #
	Guided Countersink Ø 3.25 mm	30-70325*
	Guided Countersink Ø 3.75 mm	30-70326*
	Guided Countersink Ø 4.2 mm	30-70327*
	Guided Countersink Ø 5 mm	30-70328*

Fully Guided Sleeve*

	DESCRIPTION	REF #
	Narrow, L 20 mm Drill (Brown)	30-70415
	Narrow, L 25 mm Drill (Purple)	30-70414
	Narrow, L 30 mm Drill (Silver)	30-70409
	Wide, L 20 mm Drill (Brown)	30-70413
	Wide, L 25 mm Drill (Purple)	30-70412
	Wide, L 30 mm Drill (Silver)	30-70408


Fixation Pin Sleeve*

	DESCRIPTION	REF #
	Fixation Pin Sleeve	30-70407

Fixation Pin*

	DESCRIPTION	REF #
	Fixation Pin L25mm (4 Pack)	30-70505

*Not included in the Fully Guided Surgical Kit; may be purchased separately.



Pilot Guided Drill Kit

30-70023

	DESCRIPTION	REF #
	Pilot Guided Drill Spade	30-70301
	Pilot Guided Drill Ø 2 mm L 20 mm	30-70329
	Pilot Guided Drill Ø 2 mm L 25 mm	30-70330
	Pilot Guided Drill Ø 2 mm L 30 mm	30-70331

Pilot Guided Sleeves**

	DESCRIPTION	REF #
	Ø 2 mm L 5 mm, L 20 mm Drill (Brown)	30-70405
	Ø 2 mm L 5 mm, L 25 mm Drill (Purple)	30-70403
	Ø 2 mm L 5 mm, L 30 mm Drill (Silver)	30-70401

**Not included in the Pilot Guided Drill Kit; may be purchased separately.

BASIC

KIT



Basic Complete Surgical Kit

60-70110

DESCRIPTION	QUANTITY	REF #
Pilot Drill Spade		60-70020
Step Drill Ø 1.5/2 mm, Long		60-70048
Step Drill Ø 2/2.4 mm, Long		60-70050
Step Drill Ø 2.4/3.2 mm, Long		60-70052
Step Drill Ø 3.2/3.8 mm, Long		60-70054
Final Drill for Ø 3.25 mm Implant, Long		60-70035
Final Drill for Ø 3.75 mm Implant, Long		60-70037
Final Drill for Ø 4.2 mm Implant, Long		60-70039
Final Drill for Ø 5 mm Implant, Long		60-70041
Key Ratchet Adapter	2	60-70100
Drill Extender		60-70019
Hex Key Ø 1.25 mm L 15 mm, Long		60-70101
Hex Key Ø 1.25 mm L 9 mm, Short		60-70102
Hex Key Ø 2.5 mm, Long		60-70103
Hex Key Ø 2.5 mm, Short		60-70104
Hex Key Ø 2 mm, Long		60-70105
Hex Key Ø 2 mm, Short		60-70106
Parallel Pin Ø 1.5/2 mm	4	60-70015
Manual Hex Driver Ø 1.25 mm L 14 mm, Standard		60-70006
Torque Ratchet		60-70042
Surgical Drilling Protocol		85-70041

Pilot Drill Spade

DESCRIPTION	REF #
 Pilot Drill Spade	60-70020

Step Drills, Long

DESCRIPTION	REF #
 Step Drill Ø 1.5/2 mm, Long	60-70048
 Step Drill Ø 2/2.4 mm, Long	60-70050
 Step Drill Ø 2.4/3.2 mm, Long	60-70052
 Step Drill Ø 3.2/3.8 mm, Long	60-70054

Final Drills

DESCRIPTION	REF #
 Final Drill for Ø 3.25 mm Implant, Long	60-70035
 Final Drill for Ø 3.75 mm Implant, Long	60-70037
 Final Drill for Ø 4.2 mm Implant, Long	60-70039
 Final Drill for Ø 5 mm Implant, Long	60-70041

Key Ratchet Adapter

DESCRIPTION	REF #
 Key Ratchet Adapter	60-70100

Drill Extender

DESCRIPTION	REF #
 Drill Extender	60-70019

Hex Keys, Long

DESCRIPTION	REF #
 Hex Key Ø 1.25 mm L 15 mm, Long	60-70101
 Hex Key Ø 2 mm, Long	60-70105
 Hex Key Ø 2.5 mm, Long	60-70103


Hex Keys, Short

DESCRIPTION	REF #
 Hex Key Ø 1.25 mm L 9 mm, Short	60-70102
 Hex Key Ø 2 mm, Short	60-70106
 Hex Key Ø 2.5 mm, Short	60-70104

Parallel Pins

DESCRIPTION	REF #
 Parallel Pin Ø 1.5/2 mm	60-70015

Manual Hex Driver

DESCRIPTION	REF #
 Manual Hex Driver Ø 1.25 mm L 14 mm, Standard	60-70006

Torque Ratchet


DESCRIPTION	REF #
 Torque Ratchet	60-70042



Prosthetic Kit

60-70245

	DESCRIPTION	REF #
	Hex Key Ø 1.25 mm L 6.5 mm, Extra-Short	60-70164
	Hex Key Ø 1.25 mm L 9 mm, Short	60-70102
	Hex Key Ø 1.25 mm L 15 mm, Long	60-70101
	Manual Hex Driver Ø 1.25 mm L 14 mm, Standard	60-70006
	Manual Hex Driver Ø 1.25 mm L 23 mm, Long	60-70159
	Key Ratchet Adapter	60-70100
	Torque Ratchet	60-70042



Rapid Surgical Kit

60-70149

	DESCRIPTION	QUANTITY	REF #
	Key Ratchet Adapter		60-70100
	Drill Extender		60-70019
	Hex Key Ø 1.25 mm L 15 mm, Long		60-70101
	Hex Key Ø 2.5 mm, Long		60-70103
	Hex Key Ø 2 mm, Long		60-70105
	Parallel Pin Ø 1.5/2 mm	2	60-70015
	Pilot Drill Spade		60-70020
	Step Drill Ø 2.0/2.4 mm, Long		60-70050
	Final Drill for Ø 3.25 mm Implant, Long		60-70035
	Final Drill for Ø 3.75 mm Implant, Long		60-70037
	Final Drill for Ø 4.2 mm Implant, Long		60-70039
	Final Drill for Ø 5 mm Implant, Long		60-70041
	Torque Ratchet		60-70042



Osteotome Kit

60-70080

DESCRIPTION		REF #
	Osteotome for Ø 3.25 mm Implant	60-70081
	Osteotome for Ø 3.75 mm Implant	60-70082
	Osteotome for Ø 4.2 mm Implant	60-70083
	Osteotome for Ø 5 mm Implant	60-70084
	Hand Piece Straight	60-70014



Implant Body Try In Kit

60-70142

DESCRIPTION		REF #
	Final Depth Pin for Ø 3.25 mm Implant	60-70085
	Final Depth Pin for Ø 3.75 mm Implant	60-70086
	Final Depth Pin for Ø 4.2 mm Implant	60-70087
	Final Depth Pin for Ø 5 mm Implant	60-70088



Short Step Drill Kit

60-70150

DESCRIPTION		REF #
	Step Drill Ø 1.5/2 mm, Short	60-70047
	Step Drill Ø 2.0/2.4 mm, Short	60-70049
	Step Drill Ø 2.4/3.2 mm, Short	60-70051
	Step Drill Ø 3.2/3.8 mm, Short	60-70053



Basic Osteotome Kit

60-70119

DESCRIPTION	REF #
 Osteotome for Ø 3.25 mm Implant	60-70081
 Osteotome for Ø 3.75 mm Implant	60-70082
 Osteotome for Ø 4.2 mm Implant	60-70083
 Osteotome for Ø 5 mm Implant	60-70084



Bone Tap Kit


60-70118

DESCRIPTION	REF #
 Bone Tap for Ø 3.25 mm Implant	60-70056
 Bone Tap for Ø 3.75 mm Implant	60-70057
 Bone Tap for Ø 4.2 mm Implant	60-70058
 Bone Tap for Ø 5 mm Implant	60-70059



Parallel Pin Kit

60-70115

DESCRIPTION	QUANTITY	REF #
 Final Depth Pin for Ø 3.25 mm Implant		60-70085
 Final Depth Pin for Ø 3.75 mm Implant		60-70086
 Final Depth Pin for Ø 4.2 mm Implant		60-70087
 Final Depth Pin for Ø 5 mm Implant		60-70088
 Stepped Depth Pin Ø 1.5 mm/2 mm	2	60-70134
 Stepped Depth Pin Ø 2 mm/2.4 mm	2	60-70135
 Stepped Depth Pin Ø 2.4 mm/3.2 mm	2	60-70136
 Stepped Depth Pin Ø 3.2 mm/3.8 mm	2	60-70137



Countersink Kit

60-70117

DESCRIPTION	REF #
 Countersink for Ø 5 mm & Ø 6 mm Implants	60-70143
 Countersink for Ø 4.2 mm Implant	60-70144
 Countersink for Ø 3.25 mm & Ø 3.75 mm Implants	60-70145



Bone Profiler Post Implant Kit

60-70116


DESCRIPTION	REF #
 Bone Profiler Post Implant Ø 5 mm	60-70127
 Bone Profiler Post Implant Ø 6 mm	60-70128
 Bone Profiler Post Implant Pin NP Ø 5 mm	60-70129
 Bone Profiler Post Implant Pin SP Ø 5 mm	60-70130
 Bone Profiler Post Implant Pin Conical Ø 5 mm	60-70131
 Bone Profiler Post Implant Pin SP Ø 6 mm	60-70132
 Bone Profiler Post Implant Pin Conical Ø 6 mm	60-50133



Wide Implant Drill Kit

60-70114

DESCRIPTION	REF #
 Step Drill Ø 3.8/4.3 mm, Long	60-70121
 Step Drill Ø 4.3/5.3 mm, Long	60-70123
 Countersink for Ø 5 mm & Ø 6 mm Implants	60-70143
 Final Drill for Ø 6 mm Implant, Long	60-70098



Triangular Drill Stop Kit

60-70230

DESCRIPTION	REF #
 Step Drill Stop for L 16 mm Implant, Narrow Long	60-70206
 Step Drill Stop for L 13 mm Implant, Narrow Long	60-70205
 Step Drill Stop for L 11.5 mm Implant, Narrow Long	60-70204
 Step Drill Stop for L 10 mm Implant, Narrow Long	60-70203
 Step Drill Stop for L 8 mm Implant, Narrow Long	60-70202
 Step Drill Stop for L 6 mm Implant, Narrow Long	60-70201
 Step Drill Stop for L 16 mm Implant, Wide Long	60-70212
 Step Drill Stop for L 13 mm Implant, Wide Long	60-70211
 Step Drill Stop for L 11.5 mm Implant, Wide Long	60-70210
 Step Drill Stop for L 10 mm Implant, Wide Long	60-70209
 Step Drill Stop for L 8 mm Implant, Wide Long	60-70208
 Step Drill Stop for L 6 mm Implant, Wide Long	60-70207
 Final Drill Stop for L 16 mm Implant, Extra Wide	60-70218
 Final Drill Stop for L 13 mm Implant, Extra Wide	60-70217
 Final Drill Stop for L 11.5 mm Implant, Extra Wide	60-70216
 Final Drill Stop for L 10 mm Implant, Extra Wide	60-70215
 Final Drill Stop for L 8 mm Implant, Extra Wide	60-70214
 Final Drill Stop for L 6 mm Implant, Extra Wide	60-70213


Ratchets

	DESCRIPTION	REF #
	Ratchet	60-70013
	Torque Ratchet	60-70042
	Torque Ratchet With Bending Beam	60-70079

Short Final Drills

	DESCRIPTION	REF #
	Final Drill for Ø 3 mm Implant, Short	60-70168
	Final Drill for Ø 3.25 mm Implant, Short	60-70034
	Final Drill for Ø 3.75 mm Implant, Short	60-70036
	Final Drill for Ø 4.2 mm Implant, Short	60-70038
	Final Drill for Ø 5 mm Implant, Short	60-70040
	Final Drill for Ø 6 mm Implant, Short	60-70167

Hex Drivers

	DESCRIPTION	REF #
	Manual Hex Driver Ø 1.25 mm L 9 mm, Short	60-70005
	Manual Hex Driver Ø 1.25 mm L 14 mm, Standard	60-70006
	Manual Hex Driver Ø 1.25 mm L 23 mm, Long	60-70159
	Hex Key Ø 1.25 mm L 6.5 mm, Extra-Short	60-70164
	Hex Key Ø 1.25 mm L 9 mm, Short	60-70102
	Hex Key Ø 1.25 mm L 15 mm, Long	60-70101
	Hex Key Ø 1.25 mm L 23 mm, Extra-Long	60-70165
	Hex Driver for Ratchet Ø 1.25 mm, Long	60-70003
	Hex Driver for Ratchet Ø2.5mm, Long	60-70009
	Hex Driver for Ratchet Ø2.0mm, Long	60-70012

Osteotome for DIVA

	DESCRIPTION	REF #
	Ø2.2/2.7mm	60-70090
	Ø2.7/3.2mm	60-70091
	Ø3.2/3.75mm	60-70092

Protective Abutment Apron

	DESCRIPTION	REF #
	Protective Abutment Apron	40-70130

Terms and Conditions

TERMS AND CONDITIONS OF SALE. The following are the terms and conditions by which Keystone Dental, Inc. sells its products. As used herein, the term “Keystone” means Keystone Dental, Inc. and its applicable affiliates; the term “Products” means the Keystone products and services described in this invoice, and “Purchaser” means the person or company purchasing the Products.

GOVERNING TERMS. Any purchase and shipment of Products shall be deemed to be on the terms and conditions stated herein. Any and all terms and conditions submitted by Purchaser are hereby rejected.

TAXES. The prices set forth herein do not include any sales, use, excise, ad valorem, property or other taxes applicable to the sale, use or delivery of the Products, all of which shall be paid by Purchaser separately or added to the contract price and paid by Purchaser to Keystone.

PRICES. The prices set forth for the Products are in United States Dollars. All prices shall be as of the invoice date, regardless of the date the Product is originally ordered from Keystone.

DELIVERY. Products sold hereunder shall be delivered F.O.B. shipping point. Purchaser shall assume responsibility for all subsequent delivery or shipping charges.

ACCEPTANCE. Purchaser must inspect all Products for defects and/or non-conformity within seven (7) calendar days of delivery. After this period, Purchaser will have accepted such Products.

PAYMENT. Payment terms shall be as set forth on the front of the invoice. If payment is late, the maximum statutory interest rate will apply on such late amount. Keystone may, in its sole discretion, suspend shipments of Products until full payment by Purchaser of all past due amounts to Keystone is made.

RETURN; EXCHANGE. Purchaser may return any Product within thirty (30) calendar days of the invoice date so long as there is an offsetting (exchange) order of equal or greater value. Purchaser may return any Product after thirty (30) calendar days from the invoice date, but no later than ninety (90) calendar days after the invoice date, so long as there is an offsetting (exchange) order of equal or greater value, provided, however, such returns shall be subject to a 15% re-stocking fee. Keystone will not accept any returns or exchanges after ninety (90) calendar days following the invoice date.

No biomaterial products (bone grafting and membrane) will be accepted for return or exchange. No returns or exchanges will be accepted for Products that (i) were originally sold as part of a kit, pack or set, (ii) were included on the original invoice at no charge, (iii) have less than six (6) months remaining before expiration or (iv) were custom ordered. The Purchaser shall pay all shipping costs associated with the return or exchange of a Product. No Products will be accepted by Keystone for return or exchange without an RMA number given to Purchaser by Keystone Customer Service. Returned Products must be unexpired, unopened and undamaged, returned freight prepaid.

WARRANTY. EXCEPT AS EXPRESSLY SET FORTH IN KEYSTONE’S LIFETIME IMPLANT WARRANTY PROGRAM, KEYSTONE DOES NOT MAKE ANY WARRANTY, EXPRESS OR IMPLIED, EXCEPT THAT ITS PRODUCTS SHALL BE FREE FROM DEFECTS IN MATERIAL AND/OR WORKMANSHIP. THIS WARRANTY APPLIES ONLY TO THE ORIGINAL PURCHASER. PURCHASER ASSUMES ALL RISKS AND LIABILITY RESULTING FROM THE USE OF THESE PRODUCTS, WHETHER USED SEPARATELY OR IN COMBINATION WITH OTHER PROD-

UCTS. THE WARRANTY DOES NOT APPLY IF THE FAILURE OF THE PRODUCT IS CAUSED BY TRAUMA OR ACCIDENT OR BY THE PATIENT IN WHOM THE IMPLANT WAS INSTALLED OR IF THE PURCHASER IS NOT CURRENT ON ALL OUTSTANDING PAYMENTS TO KEYSTONE AT THE TIME OF SUCH CLAIM. KEYSTONE MAY MODIFY OR TERMINATE ITS WARRANTY PROGRAM AT ANY TIME, HOWEVER ANY CHANGES TO THE WARRANTY PROGRAM WILL NOT AFFECT PRODUCTS ALREADY PURCHASED AT THE TIME OF SUCH CHANGE. PURCHASER’S SOLE REMEDY WITH RESPECT TO CLAIMS RELATED TO KEYSTONE PRODUCTS IS RESTRICTED TO THE WARRANTY POLICY.

LIMITATION OF LIABILITY. UNDER NO CIRCUMSTANCES SHALL KEYSTONE OR ANY OF ITS AFFILIATES BE LIABLE FOR BUSINESS INTERRUPTION, LOSS OF PROFITS, SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE AND FROM ANY CAUSE WHETHER BASED IN CONTRACT, WARRANTY OR TORT (INCLUDING NEGLIGENCE) OR OTHER LEGAL THEORY, EVEN IF KEYSTONE HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

INTELLECTUAL PROPERTY. All intellectual property rights in Products sold to Purchaser shall remain the exclusive property of Keystone and its licensors.

PERSONAL DATA. It is Purchaser’s responsibility to obtain patient’s consent to transfer any personal data to Keystone.

TERMINATION. Keystone may terminate this or any agreement with Purchaser to purchase its products upon immediate written notice if Purchaser fails to pay any of its debt obligations to Keystone on time.

FORCE MAJEURE. Keystone shall not be liable for any loss or damage due to failure or delay arising out of any cause beyond the reasonable control of Keystone. In the event of any failure or delay resulting from such causes, an equitable adjustment of delivery and any other appropriate terms and conditions shall be made. No such failure or delay shall be the basis for an increase in Keystone’s obligations, nor any termination by Purchaser.


VALIDITY. If any provision of these Terms and Conditions is found to be illegal or unenforceable in any respect, such illegality or unenforceability shall not affect any other provision of these Terms and Conditions, all of which shall remain enforceable in accordance with their terms.

GENERAL PROVISIONS. The purchase of the Products and these Terms and Conditions are governed by the laws of the Commonwealth of Massachusetts, without regard to conflicts of law principles. Any dispute, claim or controversy arising out of or in connection with the purchase of the Products shall be resolved by arbitration in Boston, Massachusetts. This document incorporates all oral and written representations between the parties and constitutes the entire agreement and understanding of the parties with respect to the subject matter hereof and supersedes any and all other agreements either oral or written between the parties with respect to such subject matter. No amendment or modification may be made to this document unless in writing and duly executed by an authorized representative of Keystone. Keystone may subcontract, assign or transfer its obligations or rights to a competent third party. Purchaser may not assign or transfer any of its obligations without the consent of Keystone.

IMPLANT SPECIFICATIONS




ADVANCED CLASSIC




IH	NP		SP		
	Ø 3.25 mm		Ø 3.75 mm	Ø 4.2 mm	Ø 5 mm
6 mm	—		—	20-70005	20-70011
8 mm	—		20-70017	20-70006	20-70012
10 mm	20-70018		20-70001	20-70007	20-70013
11.5 mm	20-70019		20-70002	20-70008	20-70014
13 mm	20-70020		20-70003	20-70009	20-70015
16 mm	20-70021		20-70004	20-70010	20-70016

ADVANCED+ SLIGHTLY ACTIVE APEX




IH	NP		SP			WP
	Ø 3 mm	Ø 3.25 mm	Ø 3.75 mm	Ø 4.2 mm	Ø 5 mm	Ø 6 mm
6 mm	—	—	—	20-70005P	20-70011P	20-70022P
8 mm	—	—	20-70017P	20-70006P	20-70012P	20-70023P
10 mm	20-70034P	20-70018P	20-70001P	20-70007P	20-70013P	20-70024P
11.5 mm	20-70035P	20-70019P	20-70002P	20-70008P	20-70014P	20-70025P
13 mm	20-70036P	20-70020P	20-70003P	20-70009P	20-70015P	20-70026P
16 mm	20-70037P	20-70021P	20-70004P	20-70010P	20-70016P	20-70027P

Dynamic ENHANCED STABILITY




IH	NP		SP			WP
	Ø 3 mm	Ø 3.25 mm	Ø 3.75 mm	Ø 4.2 mm	Ø 5 mm	Ø 6 mm
6 mm	—	—	—	21-70005	21-70011	21-70022
8 mm	—	—	21-70017	21-70006	21-70012	21-70023
10 mm	21-70034	21-70018	21-70001	21-70007	21-70013	21-70024
11.5 mm	21-70035	21-70019	21-70002	21-70008	21-70014	21-70025
13 mm	21-70036	21-70020	21-70003	21-70009	21-70015	21-70026
16 mm	21-70037	21-70021	21-70004	21-70010	21-70016	21-70027

PALTOP ACTIVE IMPLANT




IH	NP	SP			WP
	Ø 3.25 mm	Ø 3.75 mm	Ø 4.2 mm	Ø 5 mm	Ø 6 mm
6 mm	—	—	29-70005	29-70011	29-70022
8 mm	—	29-70017	29-70006	29-70012	29-70023
10 mm	29-70018	29-70001	29-70007	29-70013	29-70024
11.5 mm	29-70019	29-70002	29-70008	29-70014	29-70025
13 mm	29-70020	29-70003	29-70009	29-70015	29-70026
16 mm	29-70021	29-70004	29-70010	29-70016	29-70027

PALTOP ACTIVE IMPLANT TC



IH	NP	SP			WP
	Ø 3.25 mm	Ø 3.75 mm	Ø 4.2 mm	Ø 5 mm	Ø 6 mm
6 mm	—	—	29-70005TC	29-70011TC	29-70022TC
8 mm	—	29-70017TC	29-70006TC	29-70012TC	29-70023TC
10 mm	29-70018TC	29-70001TC	29-70007TC	29-70013TC	29-70024TC
11.5 mm	29-70019TC	29-70002TC	29-70008TC	29-70014TC	29-70025TC
13 mm	29-70020TC	29-70003TC	29-70009TC	29-70015TC	29-70026TC
16 mm	29-70021TC	29-70004TC	29-70010TC	29-70016TC	29-70027TC

DIVA by PALTOP



IH	SP		
	Ø 3.75 mm	Ø 4.2 mm	Ø 5 mm
11.5 mm	23-70001	23-70003	23-70005
13 mm	23-70002	23-70004	23-70006

*The implant drawing displays the dimensions for the Ø 4.2 mm x L 11.5 mm implant.

IMPLANT SPECIFICATIONS

CC Conical Connection

		CP			
		Ø 3.25 mm	Ø 3.75 mm	Ø 4.2 mm	Ø 5 mm
PCA PALTOP CONICAL ACTIVE * Ø 4.2 mm 0.5 mm 11.5 mm 7 mm Ø 1.53 mm Ø 3.75 2.1 mm	6 mm	—	—	22-70005	22-70011
	8 mm	—	22-70017	22-70006	22-70012
	10 mm	22-70018	22-70001	22-70007	22-70013
	11.5 mm	22-70019	22-70002	22-70008	22-70014
	13 mm	22-70020	22-70003	22-70009	22-70015
	16 mm	22-70021	22-70004	22-70010	22-70016
		NP		SP	
		Ø 3.25 mm	Ø 3.75 mm	Ø 4.2 mm	Ø 5 mm
Dynamic DYNAMIC * Ø 4.2 mm 0.5 mm 11.5 mm 7 mm Ø 1.4 mm Ø 3.75 2.1 mm	6 mm	—	—	24-70005	24-70011
	8 mm	—	24-70017	24-70006	24-70012
	10 mm	24-70018	24-70001	24-70007	24-70013
	11.5 mm	24-70019	24-70002	24-70008	24-70014
	13 mm	24-70020	24-70003	24-70009	24-70015
	16 mm	24-70021	24-70004	24-70010	24-70016

*The implant drawing displays the dimensions for the Ø 4.2 mm x L 11.5 mm implant.



Imperial Implant Trade Kft.

1013 Budapest Krisztina Krt. 39 B. 1. em
Telefonszám: +3670 7708507
Email: info@keystonedental.hu
www.KeystoneDental.hu